

# DNA/RNA Shield™ Fecal Collection Tube - Dx

## Instructions For Use

Cat. No. R1101-E, R1101-1-E

### Intended Use

The DNA/RNA Shield™ Fecal Collection Tube is intended for the collection and transport of clinical specimens to be analyzed by nucleic acid-based assays.

DNA/RNA Shield™ reagent pre-filled into vials ensures sample stability during transport/storage at ambient temperatures without refrigeration or specialized equipment.

The reagent is compatible with commercially available nucleic acid extraction kits and automated workflows.

### Product Description














20 x 76 mm screw cap tube prefilled with DNA/RNA Shield™ (9 mL). Spoon attached to screwcap.

### Storage & Stability

The product should be transported and stored in its original container at 5-30°C until use.

DNA/RNA Shield™ has been tested for product stability of a maximum of 24 months after manufacturing date as its shelf life.

### Label Legend

	Do not use if package is damaged		Use-by-date		Lot number
	In vitro diagnostic medical device		Catalog number		See instruction manual
	Contains sufficient for <-> tests		Manufacturer		Do-not-reuse
			Storage instruction		Unique device identifier
			Warning		European Commission Authorized Representative

### Warnings and Precautions

FOR EXTERNAL USE ONLY.

DO NOT drink, touch or remove the reagent from the vial.

The DNA/RNA Shield™ reagent can be harmful if ingested and may cause irritation if exposed to the skin and eyes.

If the contents of the tube contact your skin, wash the affected area with soap and water. If the contents of the tube are splashed in your eyes, immediately flush eyes with water. Notify your health-care provider if irritation develops. If the contents of the vial are spilled, use a new collection kit.

This product is intended for professional use or for use in POC.

DO NOT use if the product is visibly damaged.

Avoid mixing the reagent and sodium hypochlorite (bleach) or other strong acids and bases. These mixtures could release noxious gases.

Dispose of contents/container in accordance with local, regional, national or international regulations.

For specifics, consult product Safety Data Sheet (SDS) by contacting Zymo Research Corporation directly.

### Materials Required But Not Supplied (for downstream analysis)

Appropriate materials for molecular testing according to recommended protocols.



H302: Harmful if swallowed

H315: Causes skin irritation

H319: Causes eye irritation



## Instructions For Use (see diagram below)

1. Prepare and collect fecal specimen using preferred fecal specimen collection set/kit.  
Note: Method of collecting the fecal sample must prevent feces from falling into toilet water to avoid sample contamination.
2. Unscrew the collection tube cap and use the spoon to scoop **one spoonful** of feces (approximately 1 gram or 1 mL in volume) from a sample.
3. Place the sample in the collection tube.

4. Tighten the cap and shake to mix the contents thoroughly (invert 10 times) to create a suspension.  
Note: Some fecal material may be difficult to re-suspend. As long as the material is submerged, the sample is stabilized. foaming/frothing during shaking is normal.
5. Dispose of unused fecal material and thoroughly wash hands according to your institution's guidelines.



Don't let the sample go into the toilet



Collect stool into a clean container



Scoop a portion of the stool sample into the DNA/RNA Shield™ Fecal Collection Tube



Wash hands well



Zymo Research Corp.  
17062 Murphy Avenue  
Irvine, California 92614 USA

U.S. and international contact information:

Customer Support: +1 949 679 1190 ext. 1  
orders@zymoresearch.com

Technical Support: +1 949 679 1190 ext. 3  
tech@zymoresearch.com

For more contact information visit [www.zymoresearch.com](http://www.zymoresearch.com)

Copyright 2020 Zymo Research Corp. All rights reserved.



Catherine David  
42 route du périmètre  
74940 Anney, France  
+3 395 175 4280

