

ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

# OneStep qMethyl™ Panel: A method to indicate pluripotency and characterize human embryonic stem cells

Lam Nguyen, Jill Petrisko, Manuel Krispin, Xi-Yu Jia, Zymo Research Corporation, Irvine, CA


## Abstract

Pluripotency is the ability of embryonic stem cells to differentiate into multiple cell types.<sup>1</sup> Pluripotent cells have epigenetic signatures that reflect their ability to generate multiple cell types.<sup>2</sup> Different DNA methylation patterns in gene regions vary between pluripotent and differentiated cells as a result of processes such as development, carcinogenesis, genomic imprinting disorders, and cell reprogramming.<sup>3,4,5</sup> In human pluripotent cells, gene promoter regions in the **NANOG**, **RAB25**, and **PTPN6** genes have been shown to maintain low levels of DNA methylation compared to differentiated cell types.<sup>3,6</sup> Conversely, gene promoter regions of **MGMT**, **GBP3**, and **LYST** have been shown to maintain high levels of methylation in pluripotent cells compared to differentiated cell types.<sup>3,4</sup> Here we present a simple, straightforward, and bisulfite-free procedure for rapid, DNA methylation assessment for the above mentioned genes.


## Introduction

The **OneStep qMethyl™ Panel (Human Pluripotent Stem Cell Panel I)** from Zymo Research provides a simple, straightforward, and bisulfite-free procedure for rapid, DNA methylation assessment of: **NANOG**, **RAB25**, **PTPN6**, **MGMT**, **GBP3**, and **LYST** in any cell type. The reagents in the plate are already premixed and optimized for robust amplification and detection. Simply add DNA into the appropriate well and then quantitate via real-time PCR... *OneStep!*

## Workflow




The whole process is achieved in a streamlined, single, premixed reaction on a real-time, 96-well PCR plate and can be performed in < 5 hours. These specific regions of interests include **NANOG**, **RAB25**, **PTPN6**, **MGMT**, **GBP3**, and **LYST**. The performance of the panel has been validated using rigorous standards to ensure repeatability, accuracy and sensitivity. An overview of the technology is shown below:




The sample DNA is divided in two parts; a **Test Reaction** and a **Reference Reaction**. **Test Reaction** samples are digested with **Methylation Sensitive Restriction Enzymes (MSREs)** while the **Reference Reaction** samples are not (mock digest). Following digestion, DNA from both samples is used for real-time PCR.

## Assay Performance




**Figure 1. Results of OneStep qMethyl™ Panel closely correlate with actual DNA methylation percent values.** A standard curve using **NANOG**, **RAB25**, **PTPN6**, **MGMT**, **GBP3**, and **LYST** primer sets from the **OneStep qMethyl™ Panel** was performed with 0%, 25%, 50%, 75%, and 100% methylated DNA from mixtures of human non-methylated and methylated DNA. DNA samples were assayed in triplicates using the ABI™ 7500 series.



**Figure 2. Differentiated (liver cell) DNA can be distinguished from embryonic stem cell DNA using the OneStep qMethyl™ Panel.** Differentiated DNA can be detected from two DNA stem cell populations (shown in purple and orange bars) using the **OneStep qMethyl™ Panel**. Investigating the six regions of interest, liver DNA was spiked into human embryonic stem cells at 10%, 1% and 0.1%. As little as 1% (~330 cells) of liver DNA can be detected within the heterogeneous DNA population. Assay performed in duplicates on the ABI™ 7500 series.

## Biological Applications



Cell Population	RAB25	NANOG	PTPN6	MGMT	GBP3	LYST
Differentiated	+	+	+	-	-	-
Pluripotent	-	-	-	+	+	+

**Figure 3. A pluripotent stem cell line can be distinguished from a differentiated cell line using the OneStep qMethyl™ Panel.** The graph above illustrates average DNA methylation percentages for the human differentiated cells (shown as grey bars) and human stem cell line (shown as white bars) for gene regions **RAB25**, **NANOG**, **PTPN6**, **MGMT**, **GBP3**, and **LYST**. The human pluripotent stem cell line shows a distinctly different pattern than the human differentiated cells.

## Conclusions

- The **OneStep qMethyl™ Panel** can be used to efficiently quantify percent DNA methylation in regions of **RAB25**, **NANOG**, **PTPN6**, **MGMT**, **GBP3**, and **LYST**.
- The **OneStep qMethyl™ Panel** can be used to indicate pluripotency in stem cells lines by comparison with a differentiated cell line.
- The **OneStep qMethyl™ Panel** is a cost effective way to quantify DNA methylation in specific gene regions for large numbers of samples.

## References

1. Tavakoli *et al.* (2009). *BMC Cell Biology* 10 (44) doi:1186/1471-2121-10-44.
2. Meissner (2010). *Nature Biotechnology* 28 (10) doi:10.1038/nbt1684.
3. Nishino *et al.* (2011). *PLOS Genetics* 7:5 :e1002085. doi:10.1371/journal.pgen.1002085.
4. Calvanese *et al.* (2008). *PLOS One* 3:9: e3294. doi:10.1371/journal.pone.0003294.
5. Chamberlain *et al.* (2010). *PNAS* 107 (41):17668-17673.
6. Deb-Rinker *et al.* (2005). *J. Biol. Chem.* 280 (8):6257-6260.

For more information on the **OneStep qMethyl™ Panel**, **OneStep qMethyl™ Arrays**, **OneStep qMethyl™ Kit**, and **OneStep qMethyl™ Lite** go to [www.zymoresearch.com](http://www.zymoresearch.com).