

Methylated & Non-methylated pUC19 DNA Set

Standards for DNA methylation analysis workflows

Highlights

- Ideal positive and negative controls for methylation detection assays.
- *In situ* spike-in controls for monitoring bisulfite conversion efficiency in Next-Generation Sequencing experiments.
- Provided pUC19MN primer pair allows for convenient assay quality control.

Catalog Numbers:
D5017



Scan with your smart-phone camera to
view the online protocol/video.



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682

Table of Contents

Product Contents	01
Specifications.....	01
Product Description.....	02
Recommended Usage.....	03
Protocol	03
Appendix	04
Ordering Information	05
Notes.....	06
Guarantee	09

Product Contents

Methylated & Non-methylated pUC19 DNA Set	D5017	Storage Temp.
Methylated pUC19 DNA	20 ng/20 µl	-20°C
Non-methylated pUC19 DNA	20 ng/20 µl	-20°C
pUC19MN Primers	20 µl	-20°C

Specifications

Methylated pUC19 DNA

- **Source** – pUC19 plasmid purified from Dam⁻, Dcm⁻ *E. coli* [enzymatically methylated by M.SssI Methyltransferase (EC 2.1.1.37)].
- **Concentration** – 1 ng/µl in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

Non-methylated pUC19 DNA

- **Source** – pUC19 plasmid purified from Dam⁻, Dcm⁻ *E. coli*.
- **Concentration** – 1 ng/µl in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).

pUC19MN Primers

- **Concentration** – 20 µM each primer in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0).
- **Primer sequences** –

pUC19MN Forward Primer:

5' – GGTTATAGTTGTTTTTTGTGTGAAATTGTTATT – 3'

pUC19MN Reverse Primer:

5' – CTAACCTTTTACTCACATATTCTTTCCTAC – 3'

Product Description

The **Methylated & Non-methylated pUC19 DNA Set** consists of two control DNAs (methylated and non-methylated) along with a set of specifically designed primers that can be used in conjunction with the **EZ DNA Methylation™**, **EZ DNA Methylation-Gold™**, **EZ DNA Methylation-Direct™**, and **EZ DNA Methylation-Lightning™** kits from Zymo Research to assess the efficiency of bisulfite-mediated conversion of DNA. These plasmids can be used in conjunction with genomic DNA to provide internal controls to assess bisulfite conversion efficiency or to produce known mixtures of methylated and non-methylated DNA for assay calibration.

The **Non-Methylated pUC19 DNA** is pUC19 DNA that was isolated from a methylation-negative strain of *E. coli* (Dam⁻, Dcm⁻) and can be used as a negative control for DNA methylation analysis. The **Methylated pUC19 DNA** is pUC19 DNA that has been isolated from the same strain and has been enzymatically methylated at all cytosine positions comprising CG dinucleotides by M.SssI methyltransferase¹ (EC 2.1.1.37; Figure 1) and can be used as a positive control for DNA methylation analysis.




Figure 1. M.SssI methyltransferase methylates all cytosine residues in the double stranded CpG context.

The primer set herein has been designed to amplify a fragment of the supplied pUC19 DNA following bisulfite treatment. The methylated cytosines, comprising CG dinucleotides in the **Methylated pUC19 DNA** remain unconverted following bisulfite treatment, whereas nonmethylated cytosines are converted into uracil and detected as thymine after PCR. The supplied pUC19 DNA has been linearized at position 2177 using Scal endonuclease.

¹Nur et al. J. Bacteriol. 164: 19-24 (1985).

Recommended Usage

The **Methylated and Non-methylated pUC19 DNA Set** can be used in a variety of methylation analysis applications including bisulfite and methylation-specific PCR, *in situ* spike-in controls, methylation sensitive high resolution melt analysis, methylation arrays, methylated DNA immunoprecipitation (MeDIP), library preparation, and more.

Protocol

For best results, it's important to ensure the DNA is completely homogenous and fully in solution before quantification and usage. The following steps are recommended before quantification and usage¹:

1. Bring the standards to room temperature.
2. Mix the standards well, briefly spin down for 5-10 seconds.
3. Repeat Step 2 twice, for three times total.
4. Proceed with quantification or usage.

Bisulfite Conversion: For most applications 5-50 pg of plasmid may be used as an internal spike-in control for bisulfite conversion reactions containing 250 ng to 2 µg of genomic DNA. Refer to the kit specifications for setup of the bisulfite conversion reaction.

Bisulfite PCR Setup: The following is designed for a 25 µl reaction.

Component	Volume	Final Concentration
pUC19MN Primers ²	Variable	0.2 to 0.8 µM each
Bisulfite-converted DNA ³	1 µl	Up to 0.25 pg/µl
10 mM dNTP mix	0.5 µl	0.2 mM each dNTP
Standard PCR Buffer	Variable	1X
MgCl ₂ or MgSO ₄	Variable	1-4 mM, if needed
ZymoTaq™ DNA Polymerase ⁴	Variable	1-2 units
Nuclease Free Water	Bring reaction to 25 µl	N/A

Recommended Thermocycler Conditions:

- A. 95 °C, 10 minutes
- B. 95 °C, 30 seconds
- C. 57 °C, 30 to 60 seconds
- D. 72 °C, 60 seconds
- E. Repeat steps B through D an additional 29 to 39 times depending on the polymerase used.
- F. 72 °C, 7 minutes
- G. 4 °C

¹Standards are quantified using NanoDrop® measurements. If using other methods, variation may be observed.

²Alternatively, you may substitute primers of your choice.

³Remember to bisulfite-treat the DNA prior to performing PCR.

⁴We recommend using **ZymoTaq™ DNA Polymerase** or other hot-start DNA polymerases for amplification of bisulfite-treated DNA.

Appendix

pUC19MN Bisulfite PCR

The expected PCR amplicon for the Methylated & Non-methylated pUC19 DNA is 362 bp, corresponding to nucleotide positions 464 to 825 of the pUC19 sequence, including the regions (italicized) that hybridize to the primers.

Original sequence of pUC19 for bisulfite treatment and PCR amplification (sense strand 5' to 3'). The cytosines (underlined) in the CpG dinucleotide context (bold capital letters) are methylated enzymatically by M.SssI methyltransferase:

```
5' - ggtcatagct gtttcctgtg tgaaattggt atcCGctcac
aattccacac aacataCGag cCGgaagcat aaagtgtaaa gcctgggggtg
cctaatagagt gagctaactc acattaattg CGttgCGctc actgccCGct
ttccagtgCGg gaaacctgtC Gtgccagctg cattaatgaa tCGgccaaCG
CGCGgggaga ggCGgtttgC GtattgggCG ctcttcCGct tcttCGctca
ctgactCGct gCGctCGgtC GttCGgtgC GgCGagCGgt atcagctcac
tcaaaggCGg taataCGgtt atccacagaa tcaggggata aCGcaggaaa
gaacatgtga gcaaaaggcc ag - 3'
```

Expected sequence of the above PCR amplicon following bisulfite treatment:

Methylated pUC19 DNA: Below is the expected sequence for the Methylated pUC19 DNA after bisulfite conversion and PCR (sense strand). Methylated cytosines in the CpG dinucleotide context remain unconverted following bisulfite treatment, whereas non-methylated cytosines, or cytosines not in the CpG context, are converted to uracil and detected as thymine after PCR.

```
5' - ggttatagtt gttttttgtg tgaaattggt attCGtttat
aattttatat aatataCGag tCGgaagtat aaagtgtaaa gtttggggtg
tttaaatgagt gagttaattt atattaattg CGttgCGttt attggtCGtt
ttttagtgCGg gaaatttgtC Gtggttagttg tattaatgaa tCGgttaaCG
CGCGgggaga ggCGgtttgC GtattgggCG ttttttCGtt ttttCGttta
ttgattCGtt gCGttCGgtC GttCGgttgC GgCGagCGgt attagtttat
ttaaaggCGg taataCGgtt atttatagaa ttaggggata aCGtaggaaa
gaatatgtga gtaaaagggt ag - 3'
```

Non-methylated pUC19 DNA: Below is the expected sequence for the non-methylated pUC19 DNA after bisulfite conversion and PCR (sense strand). During treatment with sodium bisulfite, non-methylated cytosines are converted to uracil, which are later detected as thymine after PCR.

```
5' - ggttatagtt gttttttgtg tgaaattggt attTGtttat
aattttatat aatataTGag tTGgaagtat aaagtgtaaa gtttggggtg
tttaaatgagt gagttaattt atattaattg TGttgTGttt attggtTGtt
ttttagtgTGg gaaatttgtT Gtggttagttg tattaatgaa tTGgttaaTG
TGTGgggaga ggTGgtttgtT GtattgggTG ttttttTGtt ttttTGttta
ttgattTGtt gTGttTGgtT GttTGgttgT GgTGagTGgt attagtttat
ttaaaggTGg taataTGgtt atttatagaa ttaggggata aTGtaggaaa
gaatatgtga gtaaaagggt ag - 3'
```

Ordering Information

Product Description	Catalog No.	Size
Human Methylated & Non-methylated DNA Set	D5014	5 µg/20 µl
Human HCT116 DKO Non-Methylated DNA	D5014-1	5 µg/20 µl
Human HCT116 DKO Methylated DNA	D5014-2	5 µg/20 µl
Bisulfite-Converted Universal Methylated Human DNA Standard	D5015	1 µg/50 µl
ZymoTaq™ qPCR Premix	E2054 E2055	50 Rxns. 200 Rxns.
ZymoTaq™ Premix	E2003 E2004	50 Rxns 200 Rxns.
EZ DNA Methylation-Lightning™ Kit	D5030 D5031	50 Rxns 200 Rxns.
EZ DNA Methylation-Direct™ Kit	D5020 D5021	50 Rxns 200 Rxns.
EZ DNA Methylation™ Kit	D5001 D5002	50 Rxns 200 Rxns.
EZ DNA Methylation-Gold™ Kit	D5005 D5006	50 Rxns 200 Rxns.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are approximately 20 lines visible. The paper appears to be a standard notebook page or a sheet of stationery.



100% satisfaction guarantee on all Zymo Research products, or your money back.

Zymo Research is committed to simplifying your research with quality products and services. If you are dissatisfied with this product for any reason, please call 1(888) 882-9682.

Integrity of kit components is guaranteed for up to one year from date of purchase.
Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

™ Trademarks of Zymo Research Corporation

NanoDrop® is a registered trademark of Thermo Fisher Scientific, Inc.

The Polymerase Chain Reaction (PCR) process is covered by U.S. Patent: #4,683,195; 4,683,202 assigned to Hoffmann-La Roche. Patents pending in other countries. No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of Zymo Research's products. Further information on purchasing licenses to practice the PCR process can be obtained from the director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®

