

*For simultaneous detection of the SARS-CoV-2 virus and a host internal control*

## Instructions for Use

**For *in vitro* diagnostic (IVD) use**



### Table of Contents

|   |    |
|---|----|
| <b>Introduction</b> .....                                 | 01 |
| Intended Use .....  | 01 |
| Summary and Explanation of Test.....                      | 01 |
| Principles of the Procedure.....                          | 01 |
| <b>Reagents and Materials</b> .....                       | 02 |
| Materials Provided.....                                   | 02 |
| Reagent Storage and Handling .....                        | 02 |
| Materials Required (not provided).....                    | 03 |
| General Laboratory Warnings and Precautions .....         | 03 |
| Special Precautions.....                                  | 04 |
| Handling Precautions for Specimens .....                  | 04 |
| Work Areas .....  | 04 |
| Prevention of Nucleic Acid Contamination.....             | 04 |
| <b>Procedures</b> .....                                   | 05 |
| Quick SARS-CoV-2 Multiplex Kit Procedure.....             | 05 |
| Kit Protocol.....   | 05 |
| Amplification Reaction Set-up Procedure.....              | 05 |
| rRT-PCR Reaction Setup.....                               | 06 |
| 96-Well rRT-PCR Setup.....                                | 06 |
| 384-Well rRT-PCR Setup.....                               | 07 |
| Real-Time PCR Instrument Setup .....                      | 08 |
| Data Analysis.....  | 08 |
| Interpretation of Results.....                            | 09 |
| <b>Appendices</b> .....                                   | 10 |
| Examples of Normal and Aberrant Amplification Curves..... | 10 |
| Troubleshooting Guide.....                                | 11 |
| <b>Limitations of Procedures</b> .....                    | 12 |
| Conditions of Authorization for the Laboratory.....       | 13 |
| <b>Performance Characteristics</b> .....                  | 14 |
| <b>Symbols Legend</b> .....                               | 22 |
| <b>Technical Support</b> .....                            | 22 |



## Introduction

This package insert must be read carefully prior to use. Instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

## Intended Use

*Quick SARS-CoV-2 Multiplex Kit* is a real-time RT-PCR test intended for the qualitative detection of nucleic acids from SARS-CoV-2 in upper respiratory specimens (such as nasal, nasopharyngeal, mid-turbinate or oropharyngeal swabs), and lower respiratory specimens (such as sputum, tracheal aspirates, and bronchoalveolar lavage) from patients suspected of COVID-19 by their healthcare provider. Testing is limited to laboratories certified to perform high complexity tests.

Results are for the identification of SARS-CoV-2 RNA. The SARS-CoV-2 RNA is generally detectable in upper and lower respiratory specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2 RNA. Clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Negative results must be combined with clinical observations, patient history, and epidemiological information.

The *Quick SARS-CoV-2 Multiplex Kit* is intended for use by qualified clinical laboratory personnel specifically instructed and trained in the techniques of real-time RT-PCR assays.

## Summary and Explanation of Test

The *Quick SARS-CoV-2 Multiplex Kit* is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test on the Bio-Rad CFX96 or CFX384 Touch™ Real-Time PCR Detection System using the Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 software (or higher). The SARS-CoV-2 primer and probe set is designed to detect RNA from the SARS-CoV-2 N gene in upper and lower respiratory samples from patients who are suspected of COVID-19 by their healthcare provider.

## Principles of the Procedure

Nucleic acid from patient samples is extracted and purified as described in the procedural steps. Selective amplification of target nucleic acid from the sample is achieved by reverse transcription of the SARS-CoV-2 RNA as well as the host specific RNase P RNA and subsequent PCR amplification using the target-specific forward and reverse primers.

The *Quick SARS-CoV-2 Multiplex Kit* detects one SARS-CoV-2 specific target sequence from the Nucleocapsid gene (N) of the virus and one host specific target sequence from the RNase P gene. The primers and probes used in the rRT-PCR assay are based on the CDC published primer and probe sequences with some modification. In addition to the RNase P gene, the test uses both a no template control which is taken through all procedural steps, including the extraction, and a positive control that monitor integrity of reagents and correct performance of the testing procedure.

## Reagents and Materials

### Materials Provided

| Component Name            | Description  | Concentration  | Component Cat # | Volume | Quantity | Kit Cat # |
|---------------------------|--|--|-----------------|--------|----------|-----------|
| CV Positive Control       | SARS-CoV-2 positive control  | 2,5 copies/μl of <i>in vitro</i> transcribed SARS-CoV-2 N gene fragment spiked into 1 ng/μl human cell RNA | R3011-3-100     | 100 μl | 1        | R3013     |
|                           |  |  |                 |        | 3        | R3013-1K  |
|                           |  |  | R3011-3-3       | 3 ml   | 1        | R3013-10K |
| 2X CV Mix                 | Cocktail for one-step rRT-PCR detection of SARS-CoV-2 and RNase P. Includes: enzymes, dNTPs, MgCl <sub>2</sub> salts, additives, and the SARS-CoV-2 and RNase P primers and probes | 2X   | R3013-1-1       | 1 ml   | 1        | R3013     |
|                           |  |  |                 |        | 10       | R3013-1K  |
|                           |  |  | R3013-1-100     | 100 ml | 1        | R3013-10K |
| No Template Control (NTC) | Nuclease-free water  | N/A  | R3011-4-1       | 1 ml   | 3        | R3013     |
|                           |  |  |                 |        | 9        | R3013-1K  |
|                           |  |  | R3011-4-100     | 100 ml | 1        | R3013-10K |

Note - Integrity of kit components are guaranteed for up to the claimed expiration date on the kit under proper storage conditions. Reagents are routinely tested on a lot-to-lot basis to ensure they provide maximal performance and reliability.

### Reagent Storage and Handling

The *Quick* SARS-CoV-2 Multiplex Kit is to be shipped on dry ice.

If received in a condition other than the label indicates, or that are damaged, contact Zymo Research Corp. directly.

Upon receipt, all components of the kit should be stored at ≤ -20°C

Store at -20°C for 30 days or for longer storage, store at ≤ -70°C.

Reagents are stable for up to five (5) freeze-thaw cycles.

## Reagents and Materials Required (NOT provided):

| Product Name   | Catalog No.  | Manufacturer |
|--|--------------|--------------|
| CFX96 Touch™ Real-Time PCR Detection System (with optics capable of detecting HEX and Quasar 670) with CFX Maestro software  | 1855195      | Bio-Rad      |
| CFX384 Touch™ Real-Time PCR Detection System (with optics capable of detecting HEX and Quasar 670) with CFX Maestro software | 1855485      | Bio-Rad      |
| Microcentrifuge  | Non-specific | Non-specific |
| Mini Plate Spinner   | Non-specific | Non-specific |
| Hard Shell PCR Plate, 96-well, thin wall   | HSP9601      | Bio-Rad      |
| Hard Shell PCR Plate, 384-well, thin wall  | HSP3805      | Bio-Rad      |
| Microseal 'B' seal   | MSB1001      | Bio-Rad      |
| Aerosol barrier pipette tips (Nuclease-Free)   | Non-specific | Non-specific |
| Micropipettes (2, 10, 200, 1000 µl)  | Non-specific | Non-specific |
| Vortex mixer   | Non-specific | Non-specific |
| Freezer (≤ -70°C)  | Non-specific | Non-specific |
| Disposable gloves, powder-free   | Non-specific | Non-specific |

### General Laboratory Warnings and Precautions

This assay is for *in vitro* diagnostic use.

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect spills of specimens by including the use of soap and water (i.e., 20% aqueous solution of Sodium Dodecyl Sulfate disinfectant (SDS))
- Decontaminate and dispose of all potentially infectious materials in accordance with local, state, and federal regulations.

Important information regarding the safe handling, transport, and disposal of this product is contained in the Safety Data Sheet. Safety Data Sheets are available from Zymo Research Corp. Inquire directly.

## Special Precautions

As with any test procedure, good laboratory practice is essential to the proper performance of this test. Due to the high sensitivity of this test, care should be taken to keep reagents and amplification mixtures free of contamination.

Only personnel proficient in handling infectious materials and the use of the *Quick SARS-CoV-2 Multiplex Kit* and the Bio-Rad CFX96 or CFX384 Touch™ Real-Time PCR Detection System together with the Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 (or higher) software should perform this procedure.

### Handling Precautions for Specimens

All patient samples should be handled as if infectious, using good laboratory procedures.

Inadequate or inappropriate specimen collection, storage, and transport are likely to yield false test results. Training in specimen collection is highly recommended due to the importance of specimen quality.

Compliance with good laboratory practices is essential to minimize the risk of cross-contamination between samples and the inadvertent introduction of RNases into samples.

Amplification technologies such as PCR are sensitive to accidental introduction of product from previous amplification reactions. Incorrect results could occur if either the clinical specimen or the reagents used become contaminated by accidental introduction of even a few molecules of amplification product. Measures to reduce the risk of contamination in the laboratory include physically separating the activities involved in performing PCR in compliance with good laboratory practices.

### Work Areas

The use of 2 dedicated areas (Sample Preparation Area and Amplification Area) within the laboratory is recommended when performing the *Quick SARS-CoV-2 Multiplex Kit*.

All reagents used in the Sample Preparation Area should remain in the dedicated area at all times. Do not bring amplification product into the Sample Preparation Area.

The Amplification Area is dedicated to the amplification and detection of amplified product. Laboratory coats and equipment used in the Amplification Area must remain in this area and not be moved to the Sample Preparation Area.

Work area and instrument platforms must be considered potential sources of contamination. Change gloves after contact with potential contaminants (specimens, eluates, and/or amplified product) before handling unopened reagents, no template control, positive control, or specimens.

Decontaminate and dispose of all potentially biohazardous materials in accordance with local, national, and European regulations.

## Prevention of Nucleic Acid Contamination

Nucleic acid contamination is minimized through:

- Reverse transcription, PCR amplification, and oligonucleotide hybridization occur in a sealed 96- or 384-Well PCR Plate.
- Detection is carried out automatically without the need to open the 96-Well Reaction Plate.
- Pipettes with aerosol barrier tips are used for all pipetting, which are discarded after use.
- Separate, dedicated areas are used to perform the *Quick SARS-CoV-2 Multiplex Kit*. Refer to the **Special Precautions** section of this package insert.

## Procedure

### Quick SARS-CoV-2 Multiplex Kit Procedure

This insert contains instructions for the *Quick SARS-CoV-2 Multiplex Kit*. The *Quick SARS-CoV-2 Multiplex Kit* has been validated using upper and lower respiratory tract specimens collected in DNA/RNA Shield™ (Zymo Research, cat. R1107-E, R1109-E, R1210-E). All samples were extracted using the *Quick-DNA/RNA Viral Magbead* (Zymo Research, cat. R2140-E, R2141-E) extraction kit by both manual and automated procedures. Automated extraction was performed on the KingFisher Flex Purification System (Thermo Fisher Scientific). Real-time RT-PCR was performed using the CFX96 and CFX384 Touch Real-Time PCR Detection Systems (Bio-Rad) with CFX Maestro™ 1.1 Version 4.1.2433.1219 software.

Deviation from the validated workflow described above may result in changes to the performance of the *Quick SARS-CoV-2 Multiplex Kit*.

### Kit Protocol

Laboratory personnel should be trained to operate the Bio-Rad CFX96 or CFX384 Touch™ Real-Time PCR Detection Systems. Operators should have a thorough knowledge of the software run on the instrument/s and must follow good laboratory practices.

The Bio-Rad CFX96 or CFX384 Touch™ Real-Time PCR Detection Systems must be linked together with the Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 software and installed prior to performing the assay. For a detailed description of how to use the application refer to the software guide: <https://www.bio-rad.com/en-us/product/cfx-maestro-software-for-cfx-real-time-pcr-instruments?ID=OKZP7E15>


For a detailed description of how to use the Bio-Rad CFX96 or CFX384 Touch™ Real-Time PCR Detection System, refer to the instrument guide. (<https://www.biorad.com/webroot/web/pdf/lsr/literature/10000068706.pdf>).

### Amplification Reaction Set-up Procedure

#### Overview:

The assay can be used with RNA extracted by both manual and automated protocols.

As illustrated in the workflow figure below, extracted RNA samples are mixed directly with the 2X CV Mix and analyzed using a real-time PCR system. Results are ready to be interpreted after the reverse transcription and real-time PCR steps are complete.



**I. rRT-PCR Reaction Setup:**

**96-Well rRT-PCR Setup:**

- ✓ Before starting, thaw frozen reagents on ice, mix 10 times by inversion, centrifuge briefly, and place back on ice.
- ✓ Avoid exposing the 2X CV Mix and reactions to direct light and keep the PCR plate on ice during preparation.
- ✓ To prevent contamination, aliquot the 2X CV Mix in the 96-well PCR plate before handling any test samples or controls.

1. Add 10 µl of **2X CV Mix** to each well.
2. Add 10 µl of every RNA sample to be tested except CV Positive Control and NTC.
3. Add 10 µl of **CV Positive Control** to one dedicated well in each plate.
4. Add 10 µl of extracted **No Template Control** (NTC) to one dedicated well in each plate.

*Example plate layout for a Quick SARS-CoV-2 Multiplex Kit assay in 96-well format is illustrated below:*

|   | 1  | 2   | 3   | 4   | 5   | 6   | 7   | 8   | 9   | 10  | 11  | 12  |
|---|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | S1 | S9  | S17 | S25 | S33 | S41 | S49 | S57 | S65 | S73 | S81 | S89 |
| B | S2 | S10 | S18 | S26 | S34 | S42 | S50 | S58 | S66 | S74 | S82 | S90 |
| C | S3 | S11 | S19 | S27 | S35 | S43 | S51 | S59 | S67 | S75 | S83 | S91 |
| D | S4 | S12 | S20 | S28 | S36 | S44 | S52 | S60 | S68 | S76 | S84 | S92 |
| E | S5 | S13 | S21 | S29 | S37 | S45 | S53 | S61 | S69 | S77 | S85 | S93 |
| F | S6 | S14 | S22 | S30 | S38 | S46 | S54 | S62 | S70 | S78 | S86 | S94 |
| G | S7 | S15 | S23 | S31 | S39 | S47 | S55 | S63 | S71 | S79 | S87 | NTC |
| H | S8 | S16 | S24 | S32 | S40 | S48 | S56 | S64 | S72 | S80 | S88 | PC  |

S# = Test Sample

NTC = No Template Control

PC = CV Positive Control

5. Firmly seal the PCR plate with an optical sealing film.
6. Briefly vortex the PCR plate and centrifuge to eliminate bubbles and bring any droplets to the bottom of the well.
7. Place the PCR plate on ice.

### 384-Well rRT-PCR Setup:

- ✓ Before starting, thaw frozen reagents on ice, mix 10 times by inversion, centrifuge briefly, and place back on ice.
- ✓ Avoid exposing the 2X CV Mix and reactions to direct light and keep the PCR plate on ice during preparation.
- ✓ To prevent contamination, aliquot the 2X CV Mix in the 384-well PCR plate before handling any test samples or controls.

1. Add 5 µl of **2X CV Mix** to each well.
2. Add 5 µl of every RNA sample to be tested except CV Positive Control and NTC.
3. Add 5 µl of **CV Positive Control** to one dedicated well in each plate.
4. Add 5 µl of extracted **No Template Control** (NTC) to one dedicated well in each plate.

*Example plate layout for a Quick SARS-CoV-2 Multiplex Kit assay in 384-well format is illustrated below:*

|   | 1   | 2   | 3   | 4   | 5   | 6   | 7    | 8    | 9    | 10   | 11   | 12   | 13   | 14   | 15   | 16   | 17   | 18   | 19   | 20   | 21   | 22   | 23   | 24   |
|---|-----|-----|-----|-----|-----|-----|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|------|
| A | S1  | S17 | S33 | S49 | S65 | S81 | S97  | S113 | S129 | S145 | S161 | S177 | S193 | S209 | S225 | S241 | S257 | S273 | S289 | S305 | S321 | S337 | S353 | S369 |
| B | S2  | S18 | S34 | S50 | S66 | S82 | S98  | S114 | S130 | S146 | S162 | S178 | S194 | S210 | S226 | S242 | S258 | S274 | S290 | S306 | S322 | S338 | S354 | S370 |
| C | S3  | S19 | S35 | S51 | S67 | S83 | S99  | S115 | S131 | S147 | S163 | S179 | S195 | S211 | S227 | S243 | S259 | S275 | S291 | S307 | S323 | S339 | S355 | S371 |
| D | S4  | S20 | S36 | S52 | S68 | S84 | S100 | S116 | S132 | S148 | S164 | S180 | S196 | S212 | S228 | S244 | S260 | S276 | S292 | S308 | S324 | S340 | S356 | S372 |
| E | S5  | S21 | S37 | S53 | S69 | S85 | S101 | S117 | S133 | S149 | S165 | S181 | S197 | S213 | S229 | S245 | S261 | S277 | S293 | S309 | S325 | S341 | S357 | S373 |
| F | S6  | S22 | S38 | S54 | S70 | S86 | S102 | S118 | S134 | S150 | S166 | S182 | S198 | S214 | S230 | S246 | S262 | S278 | S294 | S310 | S326 | S342 | S358 | S374 |
| G | S7  | S23 | S39 | S55 | S71 | S87 | S103 | S119 | S135 | S151 | S167 | S183 | S199 | S215 | S231 | S247 | S263 | S279 | S295 | S311 | S327 | S343 | S359 | S375 |
| H | S8  | S24 | S40 | S56 | S72 | S88 | S104 | S120 | S136 | S152 | S168 | S184 | S200 | S216 | S232 | S248 | S264 | S280 | S296 | S312 | S328 | S344 | S360 | S376 |
| I | S9  | S25 | S41 | S57 | S73 | S89 | S105 | S121 | S137 | S153 | S169 | S185 | S201 | S217 | S233 | S249 | S265 | S281 | S297 | S313 | S329 | S345 | S361 | S377 |
| K | S10 | S26 | S42 | S58 | S74 | S90 | S106 | S122 | S138 | S154 | S170 | S186 | S202 | S218 | S234 | S250 | S266 | S282 | S298 | S314 | S330 | S346 | S362 | S378 |
| L | S11 | S27 | S43 | S59 | S75 | S91 | S107 | S123 | S139 | S155 | S171 | S187 | S203 | S219 | S235 | S251 | S267 | S283 | S299 | S315 | S331 | S347 | S363 | S379 |
| M | S12 | S28 | S44 | S60 | S76 | S92 | S108 | S124 | S140 | S156 | S172 | S188 | S204 | S220 | S236 | S252 | S268 | S284 | S300 | S316 | S332 | S348 | S364 | S380 |
| N | S13 | S29 | S45 | S61 | S77 | S93 | S109 | S125 | S141 | S157 | S173 | S189 | S205 | S221 | S237 | S253 | S269 | S285 | S301 | S317 | S333 | S349 | S365 | S381 |
| O | S14 | S30 | S46 | S62 | S78 | S94 | S110 | S126 | S142 | S158 | S174 | S190 | S206 | S222 | S238 | S254 | S270 | S286 | S302 | S318 | S334 | S350 | S366 | S382 |
| P | S15 | S31 | S47 | S63 | S79 | S95 | S111 | S127 | S143 | S159 | S175 | S191 | S207 | S223 | S239 | S255 | S271 | S287 | S303 | S319 | S335 | S351 | S367 | NTC  |
| Q | S16 | S32 | S48 | S64 | S80 | S96 | S112 | S128 | S144 | S160 | S176 | S192 | S208 | S224 | S240 | S256 | S272 | S288 | S304 | S320 | S336 | S352 | S368 | PC   |

S# = Test Sample

NTC = No Template Control

PC = CV Positive Control

5. Firmly seal the PCR plate with an optical sealing film.
6. Briefly vortex the PCR plate and centrifuge to eliminate bubbles and bring any droplets to the bottom of the well.
7. Place the PCR plate on ice.



## II. Real-Time PCR Instrument Setup

Steps of 'Real-Time PCR Machine Set Up' are to be performed in an area designated specifically for amplification reactions in order to avoid contamination of the sample processing area.

- Using the Real-Time PCR software (Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 software for the CFX96 or CFX384 Touch™ Real-Time PCR Detection System), create the following PCR protocol using the following parameters:

| Step | Temperature       | Time (min:sec) |                         |
|------|-------------------|----------------|-------------------------|
| 1    | 55 °C             | 15:00          |                         |
| 2    | 95 °C             | 10:00          |                         |
| 3    | 95 °C             | 0:05           | } 45 Cycles (Steps 3-6) |
| 4    | 72 °C             | 0:30           |                         |
| 5    | 57 °C             | 0:30           |                         |
| 6    | <b>Plate Read</b> |                |                         |
| 7    | <b>END</b>        |                |                         |

- Assign each well the corresponding sample names, targets, and fluorophores. 2X CV Mix contains two different fluorophores.

| Target     | Fluorophores | Ex.    | Em.    |
|------------|--------------|--------|--------|
| SARS-CoV-2 | HEX™         | 535 nm | 556 nm |
| RNase P    | Quasar® 670  | 647 nm | 670 nm |

- Enable all filters required to detect both fluorophores.
- Load the PCR plate into the CFX96 or CFX384 Touch™ Real-Time PCR Detection System (Bio-Rad) and start the PCR protocol

## III. Data Analysis (Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 software)

The Data Analysis steps described below are for the CFX96 and CFX384 Touch™ Real-Time PCR Detection System (Bio-Rad) with the CFX Maestro Software 1.1 Version 4.1.2433.1219. Adjustments to the procedure may be necessary to perform the reaction using a different real-time PCR instrument and software.

- Visually inspect the plate for any issues that occurred during PCR (e.g. evaporation) and take note of any problems.
- Adjust the Baseline Threshold value to 50 RFU for each fluorophore. Under "Settings", click on "Baseline Threshold", select "User Defined" and enter 50 into the user defined field.
- Under "Settings", click on "Baseline Setting", and select "Apply Fluorescence Drift Correction."
- Note and record any samples, PC, or NTC with amplification in the HEX™ or Quasar® 670 channel before 10 cycles. This step is necessary to identify samples with Ct values < 10 which may be excluded by adjustments made in Step 5.
- Exclude the first 10 cycles of the PCR by going to "Settings", clicking on "Cycles to Analyze". Manually enter the range as 10 to 45.

6. If aberrant amplification curves are observed, please refer to the Appendix on page 10 for examples and detailed explanations.

#### IV. Interpretation of Results

Guidelines for the Interpretation of Controls written in this section are based upon results obtained using a CFX96 or CFX384 Touch™ Real-Time PCR Detection System (Bio-Rad). Ct values provided in this section are representative and may vary between different real-time PCR instruments.

All controls should be examined prior to the interpretation of the test samples:

**CV Positive Control (PC):** The CV Positive Control will show amplification signals in both the HEX™ and Quasar® 670 channels, which detect the presence of SARS-CoV-2 and RNase P, respectively. Signals generated from the CV Positive Control are considered valid if the cycle threshold (Ct) value is  $\leq 40$  for the virus target (HEX™ channel) and  $\leq 30$  for RNase P (Quasar® 670 channel) when using the recommended systems settings. If signals for the CV Positive Control are detected after 40 amplification cycles ( $Ct > 40$ ) for the virus target or after 30 amplification cycles ( $Ct > 30$ ) for RNase P, the control must be replaced with a new aliquot. If this problem is not resolved, the whole kit must be replaced with a new one.


**No Template Control (NTC):** The No Template Control sample will show no amplification signal for either the virus target (HEX™ channel) or RNase P (Quasar 670® channel). The No Template Control sample will be considered valid if no amplification occurs in either the HEX™ or Quasar® 670 channels. If amplification occurs in either channel, contamination of extraction and/or rRT-PCR reagents may have occurred, and reagents must be replaced.

| Sample              | SARS-CoV-2 (HEX™)             | RNase P (Quasar® 670)            | Interpretation |
|---------------------|-------------------------------|----------------------------------|----------------|
| No Template Control | No amplification              | No amplification                 | Valid          |
| CV Positive Control | $Ct \leq 40$                  | $Ct \leq 30$                     | Valid          |
| Test Sample         | $Ct \leq 40$                  | Any amplification                | Positive       |
|                     | No amplification or $Ct > 40$ | $Ct < 40$                        | Negative       |
|                     | No amplification              | No amplification or $Ct \geq 40$ | Invalid        |


## Appendices

### Examples of Normal and Aberrant Amplification Curves


Note that the first 10 cycles of PCR have been excluded from the graphs, the baseline threshold for both fluorophores have been adjusted to 50 RFU, and fluorescence drift correction has been applied as described on page 8-9.




**Figure 1.** Curves A through E are examples of normal amplification curves for the SARS-CoV-2 target. Curve F (red) represents the SARS-CoV-2 signal in the CV Positive Control. Curve G (blue) represents the No Template Control.



**Figure 2.** Examples of samples with normal amplification curves for the IC (RNase P) are shown in brown. The IC curve for the CV Positive Control is indicated in red (A) and the No Template Control in blue (B).



**Figure 3.** Curves A, B, and C are examples of weak signals for the SARS-CoV-2 target that would not be considered a positive test result ( $C_t > 40$ ). For these samples, a repeat rRT-PCR may be performed to confirm the results. Note that the Y axis scale has been adjusted (magnified) to better visualize the difference between curves.



**Figure 4.** Curve A is an example of an aberrant amplification for the SARS-CoV-2 target. Curves B and C are examples of aberrant or weak positive amplification. In cases like these, results should be confirmed by repeat testing. The Y axis scale has been adjusted (magnified) to better visualize details of the curves.

## Troubleshooting Guide

| Problem  | Possible Causes and Suggested Actions  |
|--|--|
| <p><i>Reduction in volume observed in the wells after rRT-PCR.</i></p>   | <ul style="list-style-type: none"> <li>• <b>Cause:</b> Evaporation of rRT-PCR reaction mixture due to improper sealing of the PCR plate.</li> <li>• <b>Action:</b> Record the wells with reduced volume and proceed with normal analysis. Repeat rRT-PCR on samples experiencing evaporation if results are invalid.</li> </ul>  |
| <p><i>Aberrant amplification after pre-analytical steps have been executed (page 8-9). Please see examples of aberrant and normal amplification curves on page 10.</i></p> | <ul style="list-style-type: none"> <li>• <b>Cause:</b> Presence of air bubbles in the reaction, poor quality sample, incorrect setup, or compromised reagents.</li> <li>• <b>Action:</b> Record the wells and fluorophores that show aberrant amplification. Samples that exhibit aberrant amplification should be repeated to determine a conclusive result. Re-extraction of the sample or use of new reagents may be necessary if re-testing the sample does not produce a clear signal.</li> </ul>             |
| <p><i>SARS-CoV-2 signal in CV Positive Control detected after 40 Cycles.</i></p>   | <ul style="list-style-type: none"> <li>• <b>Cause:</b> Incorrect rRT-PCR set-up or the <i>Quick SARS-CoV-2 Multiplex Kit</i> reagents may have been compromised (e.g. improper storage or more than 5 freeze-thaw cycles).</li> <li>• <b>Action:</b> Replace the control. If the problem persists, the entire kit should be replaced.</li> </ul>   |
| <p><i>RNase P signal in CV Positive Control detected after 30 Cycles.</i></p>  | <ul style="list-style-type: none"> <li>• <b>Cause:</b> Incorrect rRT-PCR setup or the <i>Quick SARS-CoV-2 Multiplex Kit</i> reagents may have been compromised (e.g. improper storage or more than 5 freeze-thaw cycles).</li> <li>• <b>Action:</b> Replace the control. If the problem persists, the entire kit should be replaced.</li> </ul>  |
| <p><i>SARS-CoV-2 and/or RNase P amplification in the No Template Control (NTC)</i></p>   | <ul style="list-style-type: none"> <li>• <b>Cause:</b> Contamination from the environment, contamination of extraction and/or rRT-PCR reagents, or well-to-well cross contamination.</li> <li>• <b>Action:</b> Clean and decontaminate all surfaces and instruments. Ensure that filter tips are used during the procedure and changed between samples. Ensure that extraction and rRT-PCR setup is properly executed. Repeat extraction and rRT-PCR with new reagents including all required controls.</li> </ul> |

## Limitations of Procedures

Testing is limited to laboratories certified to perform high complexity tests.

Optimum performance of this kit requires upper and lower respiratory tract specimens collected in DNA/RNA Shield™ (Zymo Research, Cat. R1107, R1109, R1124, R1210), extracted using the *Quick*-DNA/RNA Viral MagBead Kit (Zymo Research, Cat. R2140, R2141), both manually or automated (using the KingFisher Flex; Thermo Fisher Scientific), and analyzed using the Bio-Rad CFX96 Touch™ and the CFX384 Touch™ Real-Time PCR Detection Systems and the Bio-Rad CFX Maestro™ 1.1 Version 4.1.2433.1219 software (Bio-Rad). Variations in any component of the validation workflow indicated above may result in changes in performance characteristics of the *Quick* SARS-CoV-2 Multiplex Kit.

Detection of SARS-CoV-2 RNA may be affected by sample collection methods, patient factors (e.g. presence of symptoms), and/or stage of infection.

False-negative results may arise from improper collection, shipping, and/ or storage of specimen.

Pooled sample testing may decrease the sensitivity of the *Quick* SARS-CoV-2 Multiplex Kit since it lowers the volume of each individual sample that is analyzed.

Performance of the test was established in sputum. A shift in Ct was observed for all targets at high mucin concentrations. Therefore, high mucin concentrations at or above 0.1% ( $\geq 1$  mg/ml) may result in invalid results.

The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutic or immunosuppressant drugs have not been evaluated.

As with any molecular test, mutations within the target regions detected by the *Quick* SARS-CoV-2 Multiplex Kit could affect primer and/or probe binding resulting in failure to detect the presence of virus.

The performance of this SARS-CoV-2 assay was established using sputum. Nasal, nasopharyngeal, oropharyngeal and mid-turbinate swabs, BAL and tracheal aspirates are also considered acceptable specimen types for use with the SARS-CoV-2 Assay but performance has not been established. Testing of nasal and mid-turbinate nasal swabs (self-collected at a healthcare site or collected by a healthcare provider) is limited to patients with symptoms of COVID-19. Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.

Negative results do not preclude infection with the SARS-CoV-2 virus and should not be the sole basis of a patient treatment/management or public health decision. Follow up testing should be performed according to the current Health authorities' recommendations.

- Members of the Infectious disease laboratory will be trained to perform this assay and competency will be assessed and documented.
- A false negative result may occur if a specimen is improperly collected, transported, or handled. False negative results may also occur if amplification inhibitors are present in the specimen or if inadequate numbers of organisms are present in the specimen.
- This test cannot rule out diseases caused by other bacterial or viral pathogens.

**Conditions of Authorization for the Laboratory**

- A. Laboratories using your product will include with result reports of your product, all authorized Fact Sheets. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- B. Laboratories using your product will have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- C. Laboratories will collect information on the performance of your product and report to Health Authorities and Zymo Research any suspected occurrence of false positive or false negative results and significant deviations from the established performance characteristics of your product of which they become aware.
- D. All laboratory personnel using your product must be appropriately trained in RT-PCR techniques and use appropriate laboratory and personal protective equipment when handling this kit and use your product in accordance with the authorized labeling.

## Performance Characteristics

This insert contains instructions for the *Quick SARS-CoV-2 Multiplex Kit*. The *Quick SARS-CoV-2 Multiplex Kit* has been validated using upper and lower respiratory tract specimens collected in DNA/RNA Shield™ (Zymo Research, cat. R1107-E, R1109-E, R1210-E). All samples were extracted using the *Quick-DNA/RNA Viral Magbead* (Zymo Research, cat. R2140-E, R2141-E) extraction kit by both manual and automated procedures. Automated extraction was performed on the KingFisher Flex Purification System (Thermo Fisher Scientific). Real-time RT-PCR was performed using the CFX96 and CFX384 Touch Real-Time PCR Detection System (Bio-Rad) with CFX Maestro™ 1.1 Version 4.1.2433.1219 software.

Deviation from the validated workflow described above may result in changes to the performance of the *Quick SARS-CoV-2 Multiplex Kit*.

Detailed information about the validation workflow and protocols are available upon request at [tech@zymoresearch.com](mailto:tech@zymoresearch.com)

### I. Limit of Detection (LoD) – Analytical Sensitivity

The LoD of the *Quick SARS-CoV-2 Multiplex Kit* utilizing the entire test system from sample preparation to detection, has been determined. A preliminary LoD, was determined using heat-inactivated SARS-CoV-2, strain USA-WA1/2020 (VR-1986HK, ATCC) spiked into sputum collected using the DNA/RNA Shield™ Saliva/Sputum Collection Kit. Heat-inactivated virus was quantified by ATCC using quantitative RT-PCR and reported as copies per milliliter (copies/ml). The preliminary LoD study was performed using automated nucleic acid extraction of sputum samples containing 10-fold dilutions of virus from  $8.33 \times 10^4$  to 0 copies/ml. All samples were processed using the *Quick-DNA/RNA Viral Magbead* (Zymo Research Corp.) extraction kit automated on the KingFisher Flex Purification System (Thermo Fisher Scientific). Real-time RT-PCR was performed using the Bio-Rad CFX96 Touch Real-Time PCR Detection System. The preliminary LoD was determined to be 83 copies/ml (5 copies/rxn), the lowest concentration for which 5/5 independent replicates tested positive for SARS-CoV-2 (**Table 1**).

**Table 1: Preliminary LoD study in sputum specimens using automated extraction**

| Concentration in Dilution Tested [copies/ml] | Concentration in Dilution Tested [copies/ <u>rxn</u> ] | Replicate 1 C <sub>t</sub> | Replicate 2 C <sub>t</sub> | Replicate 3 C <sub>t</sub> | Replicate 4 C <sub>t</sub> | Replicate 5 C <sub>t</sub> | Call Rate | Average Ct | Lowest Concentration with Uniform Positivity | Preliminary Limit of Detection (LoD)     |
|--|--|----------------------------|----------------------------|----------------------------|----------------------------|----------------------------|-----------|------------|--|--|
| $8.33 \times 10^4$                           | 5,000  | 25.55                      | 25.68                      | 25.63                      | 25.82                      | 25.59                      | 5/5       | 25.68      | 83 copies/ml<br>(5 copies/ <u>rxn</u> )      | <u>83 GEC/ml</u><br>(5 GEC/ <u>rxn</u> ) |
| $8.33 \times 10^3$                           | 500  | 29.11                      | 28.97                      | 29.32                      | 28.87                      | 29.10                      | 5/5       | 29.06      |  |  |
| $8.33 \times 10^2$                           | 50   | 31.87                      | 32.57                      | 32.05                      | 32.12                      | 32.88                      | 5/5       | 32.40      |  |  |
| $8.3 \times 10^1$                            | 5  | 37.07                      | 36.22                      | 35.58                      | 35.05                      | 37.71                      | 5/5       | 36.14      |  |  |
| $8.3 \times 10^0$                            | 0.5  | NA*                        | NA                         | 39.35                      | 38.65                      | NA                         | 2/5       | 39.00      |  |  |

\*NA: No amplification

The LoD of the assay was confirmed for both manual and automated nucleic acid extraction methods using sputa as a clinical matrix. Heat-inactivated SARS-CoV-2 was spiked into negative sputa samples and 20 replicates were independently processed for LoD determination. The lowest concentration for which all 5 replicates were positive in the preliminary LoD evaluation (i.e. 83 copies/ml) was used as a starting point for the confirmatory LoD study. Testing of 83 copies/mL could not confirm the tentative LoD; therefore, a concentration 2-fold above 83 copies/mL was also tested. The final LoD for which  $\geq 19/20$  replicates tested positive for both manual and automated nucleic acid extraction methods was determined to be 167 copies/ml (10 copies/rxn) (**Table 2**). Results are summarized in **Table 3**.

**Table 2: Confirmatory LoD for sputum specimens**

| Replicate | Automated Extraction          |                             |                               |   |                                       | Manual Extraction |                               |   |                                       |
|-----------|-------------------------------|-----------------------------|-------------------------------|---|---------------------------------------|-------------------|-------------------------------|---|---------------------------------------|
|           | 167 copies/ml (10 copies/rxn) | 83 copies/ml (5 copies/rxn) | 42 copies/ml (2.5 copies/rxn) | Lowest Concentration with at least 19/20 (95%) Positive | Confirmatory Limit of Detection (LoD) | Replicate         | 167 copies/ml (10 copies/rxn) | Lowest Concentration with at least 19/20 (95%) Positive | Confirmatory Limit of Detection (LoD) |
|           | Ct                            | Ct                          | Ct                            |   |                                       |                   | Ct                            |   |                                       |
| 1         | 34.30                         | 40.24                       | 38.71                         | 167 copies/ml (10 copies/rxn)                           | 167 copies/ml (10 copies/rxn)         | 1                 | 35.05                         | 167 copies/ml (10 copies/rxn)                           | 167 copies/ml (10 copies/rxn)         |
| 2         | 35.24                         | 39.55                       | 36.09                         |   |                                       | 2                 | 36.22                         |   |                                       |
| 3         | 36.05                         | 36.02                       | 37.10                         |   |                                       | 3                 | 35.04                         |   |                                       |
| 4         | 36.24                         | 36.67                       | 38.25                         |   |                                       | 4                 | 33.75                         |   |                                       |
| 5         | 35.57                         | 38.22                       | 40.26                         |   |                                       | 5                 | 35.82                         |   |                                       |
| 6         | 35.10                         | 36.58                       | 37.83                         |   |                                       | 6                 | 35.57                         |   |                                       |
| 7         | 33.85                         | 35.01                       | 36.92                         |   |                                       | 7                 | 35.40                         |   |                                       |
| 8         | 37.78                         | 37.71                       | 37.71                         |   |                                       | 8                 | 34.76                         |   |                                       |
| 9         | 35.20                         | NA                          | 36.69                         |   |                                       | 9                 | 36.08                         |   |                                       |
| 10        | 36.59                         | 36.26                       | NA                            |   |                                       | 10                | 35.73                         |   |                                       |
| 11        | 34.58                         | 35.18                       | 37.93                         |   |                                       | 11                | 38.07                         |   |                                       |
| 12        | 33.41                         | 37.72                       | 37.95                         |   |                                       | 12                | 35.75                         |   |                                       |
| 13        | 34.59                         | 35.28                       | 37.45                         |   |                                       | 13                | 34.85                         |   |                                       |
| 14        | 34.61                         | 35.92                       | NA                            |   |                                       | 14                | 35.08                         |   |                                       |
| 15        | 34.45                         | 35.45                       | 37.92                         |   |                                       | 15                | 35.24                         |   |                                       |
| 16        | 35.37                         | 35.27                       | 36.98                         |   |                                       | 16                | 35.77                         |   |                                       |
| 17        | 35.34                         | 35.33                       | 37.65                         |   |                                       | 17                | 35.60                         |   |                                       |
| 18        | 34.35                         | 39.84                       | 39.50                         |   |                                       | 18                | 34.99                         |   |                                       |
| 19        | 34.55                         | 35.68                       | 36.12                         |   |                                       | 19                | 34.20                         |   |                                       |
| 20        | 33.95                         | NA                          | 38.94                         |   |                                       | 20                | 34.48                         |   |                                       |
| Call Rate | 20/20                         | 17/20                       | 17/20                         |   |                                       | Call Rate         | 20/20                         |   |                                       |

**Table 3: Confirmatory LoD Study - Summary**

| Target Level                  | Extraction Method | Valid results | SARS-CoV-2 |            |                | Internal Control |            |                |
|-------------------------------|-------------------|---------------|------------|------------|----------------|------------------|------------|----------------|
|                               |                   |               | Positive   |            |                | Positive         |            |                |
|                               |                   |               | n          | Average Ct | Detection Rate | n                | Average Ct | Detection Rate |
| 167 copies/ml (10 copies/rxn) | Automated         | 20            | 20         | 35.06      | 100%           | 20               | 23.47      | 100%           |
| 83 copies/ml (5 copies/rxn)   | Automated         | 20            | 17         | 36.77      | 85%            | 20               | 23.35      | 100%           |
| 42 copies/ml (2.5 copies/rxn) | Automated         | 20            | 17         | 37.78      | 85%            | 20               | 23.41      | 100%           |
| 167 copies/ml (10 copies/rxn) | Manual            | 20            | 20         | 35.37      | 100%           | 20               | 23.18      | 100%           |



## II. Clinical Evaluation

A clinical evaluation was performed using 30 SARS-CoV-2 positive and 30 negative sputum specimens collected using the DNA/RNA Shield™ Saliva/Sputum Collection Kit. Specimens were collected from patients experiencing one or more symptoms of respiratory infection under the supervision of a healthcare provider. SARS-CoV-2 status was determined using the *Quick* SARS-CoV-2 Multiplex Kit and a comparator rRT-PCR assay authorized by the FDA for emergency use. Results are summarized in **Table 4**.

**Table 4: Clinical Evaluation Summary**

| Comparator Result | Quick SARS-CoV-2 Multiplex Kit |          | Total | Percent Agreement    |
|-------------------|--------------------------------|----------|-------|----------------------|
|                   | Positive                       | Negative |       | 95% CI               |
| Positive          | 30                             | 0        | 30    | PPA: 100% (30/30)    |
|                   |                                |          |       | 95% CI: 88.7% - 100% |
| Negative          | 0                              | 30       | 30    | NPA: 100% (30/30)    |
|                   |                                |          |       | 95% CI: 88.7% - 100% |
| Total             | 30                             | 30       | 60    |                      |

## III. Transport Medium Study

Additional transport medias commonly used for SARS-CoV-2 testing were evaluated for their compatibility for use with *Quick* SARS-CoV-2 Multiplex Kit compared to DNA/RNA Shield™. Performance was evaluated using sputum spiked into each transport medium (listed in **Table 5**) to a final concentration of 20%, which represents a complex clinical matrix equivalent to a swab specimen. SARS-CoV-2 low positive samples were created by spiking in 334 copies/ml (2X LoD) of heat-inactivated SARS-CoV-2 into the clinical matrix. All samples were extracted using an automated method. Ct values for both SARS-CoV-2 and the internal control for each transport medium are shown in **Table 6**. The average Ct values for each transport medium were within 1 Ct value of DNA/RNA Shield™, indicating that all transport medias tested are compatible for use with the *Quick* SARS-CoV-2 Multiplex Kit.

**Table 5. List of transport medias tested**

| Transport Medium             | Manufacturer   |
|------------------------------|--|
| DNA/RNA Shield™              | Zymo Research  |
| Viral Transport Medium (VTM) | Prepared in-house according CDC instructions ( <a href="https://www.cdc.gov/coronavirus/2019-ncov/downloads/Viral-Transport-Medium.pdf">https://www.cdc.gov/coronavirus/2019-ncov/downloads/Viral-Transport-Medium.pdf</a> ) |
| Saline                       | Prepared in-house (0.9% NaCl in nuclease-free water)   |
| PBS                          | Gibco  |
| Liquid Amies                 | Copan  |
| PrimeStore MTM               | Longhorn   |

**Table 6. Results from SARS-CoV-2 negative and low positive samples collected in different transport medias**

| SARS-CoV-2                | Transport Media | SARS-CoV-2     |                |                |                |                |            |               | Internal Control |                |                |                |                |            |               |
|---------------------------|-----------------|----------------|----------------|----------------|----------------|----------------|------------|---------------|------------------|----------------|----------------|----------------|----------------|------------|---------------|
|                           |                 | Replicate 1 Ct | Replicate 2 Ct | Replicate 3 Ct | Replicate 4 Ct | Replicate 5 Ct | Average Ct | Std Deviation | Replicate 1 Ct   | Replicate 2 Ct | Replicate 3 Ct | Replicate 4 Ct | Replicate 5 Ct | Average Ct | Std Deviation |
| Negative                  | DNA/RNA Shield  | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 25.35            | 25.24          | 25.31          | 25.26          | 25.41          | 25.31      | 0.06          |
|                           | VTM             | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 27.15            | 27.33          | 27.18          | 27.10          | 27.18          | 27.19      | 0.07          |
|                           | Saline          | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 25.00            | 24.67          | 24.68          | 24.95          | 24.72          | 24.80      | 0.14          |
|                           | PBS             | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 24.80            | 24.73          | 24.70          | 24.69          | 24.86          | 24.75      | 0.06          |
|                           | Liquid Amies    | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 25.35            | 25.39          | 25.16          | 25.44          | 25.42          | 25.35      | 0.10          |
|                           | Probes/MTM      | NA             | NA             | NA             | NA             | NA             | NA         | NA            | 26.52            | 26.47          | 26.57          | 26.56          | 26.69          | 26.56      | 0.07          |
| 2X $10^4$ (334 copies/ml) | DNA/RNA Shield  | 33.32          | 33.72          | 33.27          | 34.00          | 33.36          | 33.53      | 0.28          | 25.56            | 25.51          | 25.36          | 25.23          | 25.29          | 25.39      | 0.13          |
|                           | VTM             | 32.81          | 33.15          | 34.23          | 32.90          | 33.29          | 33.28      | 0.51          | 27.14            | 26.80          | 27.00          | 27.01          | 26.82          | 26.95      | 0.13          |
|                           | Saline          | 34.15          | 33.83          | 34.16          | 33.17          | 33.08          | 33.68      | 0.47          | 24.61            | 24.51          | 24.55          | 24.47          | 24.44          | 24.52      | 0.06          |
|                           | PBS             | 33.67          | 33.45          | 33.52          | 33.69          | 33.63          | 33.59      | 0.09          | 25.22            | 24.77          | 24.77          | 24.81          | 25.09          | 24.93      | 0.19          |
|                           | Liquid Amies    | 33.95          | 33.84          | 33.54          | 33.33          | 34.02          | 33.74      | 0.26          | 25.33            | 25.41          | 25.27          | 25.12          | 25.14          | 25.25      | 0.11          |
|                           | Probes/MTM      | 34.03          | 33.17          | 33.46          | 33.49          | 33.57          | 33.54      | 0.28          | 26.73            | 26.56          | 26.56          | 26.50          | 26.35          | 26.54      | 0.12          |

**VI. Analytical Sensitivity and Specificity**

The studies below were performed during validation of the *Quick SARS-CoV-2 rRT-PCR Kit* (Zymo Research, cat. R3011), which contains components and primer and probe sequences identical to the *Quick SARS-CoV-2 Multiplex Kit* (cat. R3013). The primers and probes used in both kits are listed above in **Table 7**. Data shown in **Tables 8 – 16** include results from only the primers and probe which are included in the *Quick SARS-CoV-2 Multiplex Kit*.

**Inclusivity (Analytical Sensitivity)**

In silico inclusivity analysis of the oligo sets for SARS-CoV-2 (taxonomy ID 2697049) was performed using multiple sequence alignment of the individual SARS-CoV-2 primers and probe against all SARS-CoV-2 N gene sequences found in the NCBI database (1,354 in total as of April 24, 2020). The primers and probe correspond to CDC’s primers and probe with unmodified sequence. Each primer and probe was found to have single nucleotide mismatches within a small number ( $\leq 0.74\%$ ) of the N gene sequences analyzed (**Table 8**). None of the mismatches occurred within the last five 3’ nucleotides of any of the primers, reducing the potential detrimental effect of those mismatches on PCR efficiency. Overall, the number of sequences with unfavorable mismatches is too small to impact the inclusivity (analytical sensitivity) of the test.

**Table 7. Primers and probes used in the *Quick SARS-CoV-2 rRT-PCR Kit* and *Quick SARS-CoV-2 Multiplex Kit*. All primers and probes with the same name have identical sequence.**

| Quick SARS-CoV-2 rRT-PCR Kit<br>(cat. R3011) | Quick SARS-CoV-2 Multiplex Kit<br>(cat. R3013) |
|--|--|
| <b>CV Mix 1 (SARS-CoV-2)</b>                 | <b>2X CV Mix (SARS-CoV-2 and IC)</b>           |
| <b>SARS-CoV-2 Target 1</b>                   | <b>SARS-CoV-2 Target</b>                       |
| CV1_Fw                                       | CV1_Fw   |
| CV1_Rv                                       | CV1_Rv   |
| CV1_PB (HEX)                                 | CV1_PB (HEX)                                   |
| <b>SARS-CoV-2 Target 2</b>                   | <b>IC Target (RNase P)</b>                     |
| CV2_Fw                                       | RP-Fw  |
| CV2_Rv                                       | RP_Rv  |
| CV2_PB (HEX)                                 | RP_PB (Quasar 670)                             |
| <b>SARS-CoV-2 Target 3</b>                   |  |
| CV3_Fw                                       |  |
| CV3_Rv                                       |  |
| CV3_PB (HEX)                                 |  |
| <b>CV Mix 2 (IC, RNase P)</b>                |  |
| RP-Fw  |  |
| RP_Rv  |  |
| RP_PB (Quasar 670)                           |  |

**Table 8. Primer and probe single nucleotide mismatches to SARS-CoV-2 N gene sequences.** SARS-CoV-2 N gene accession numbers in bold have single nucleotide mismatches that occur at the same position within the primer/probe.

| Primer/Probe | Total Number of SARS-CoV-2 N Gene Sequences Analyzed | Number of N Gene Sequences with 1 Mismatch | % N Gene Sequences with 1 Mismatch | Number of N Gene Sequences Having > 1 Mismatch | NCBI Accession Numbers for N Genes with 1 Mismatch  |
|--------------|--|--|------------------------------------|--|---|
| CV1_Fw       | 1,354  | 2  | 0.15                               | 0  | MT293178, MT350243  |
| CV1_Rv       | 1,354  | 1  | 0.07                               | 0  | MT293178  |
| CV1_PB       | 1,354  | 10   | 0.74                               | 0  | MT371038, MT326026, <b>MT372482</b> , <b>MT372481</b> , <b>MT372480</b> , <b>MT344946</b> , <b>MT304476</b> , <b>MT304475</b> , <b>MT304474</b> , <b>MT293161</b> |

### Cross-reactivity (Analytical Specificity)

Cross reactivity studies were performed using the NCBI Basic Alignment Search Tool (BLAST) to identify the largest regions of homology between the primers and probe to the genomes indicated (**Table 9**). Homology  $\geq 80\%$  was found for one primer with *Neisseria meningitidis*. A microbial interference study using samples containing  $1 \times 10^6$  *Neisseria meningitidis* cells with and without 3X LoD (750 GEC/ml) of whole genome SARS-CoV-2 RNA showed no interference (false positive or false negative result) from this microorganism with the test (**Table 10**).

*In silico* analysis showed homology of one primer and probe to the closely related SARS-CoV; however, since one primer has only 40% homology to SARS-CoV, the primers and probe used in this test are not expected to detect SARS-CoV. The primers and probe used in this test were found to have little cross-reactivity with endemic coronaviruses (229E, NL63, HKU1, and OC43) or MERS-CoV.

**Table 9. *In silico* cross-reactivity analysis of primers and probe**

| Pathogen                        | Strain                      | GenBank Accession # | % Homology |        |        |
|---------------------------------|-----------------------------|---------------------|------------|--------|--------|
|                                 |                             |                     | CV1_FW     | CV1_RV | CV1_PB |
| Human coronavirus               | 229E                        | NC_002645.1         | 50         | 46     | 38     |
|                                 | OC43                        | NC_006213.1         | 50         | 38     | 42     |
|                                 | HKU1                        | NC_006577.2         | 50         | 42     | 42     |
|                                 | NL63                        | NC_005831.2         | 45         | 42     | 42     |
| SARS coronavirus                | SARS                        | NC_004718.3         | 40         | 92     | 92     |
| MERS coronavirus                | MERS                        | NC_019843.3         | 50         | 42     | 38     |
| Human adenovirus                | Human adenovirus 1          | AC_000017.1         | 45         | 38     | 38     |
| Human metapneumovirus           | 00-1                        | NC_039199.1         | 45         | 42     | 33     |
| Human <del>respirovirus 1</del> | Human parainfluenza virus 1 | NC_003461.1         | 55         | 38     | 38     |
| Human <del>rubulavirus 2</del>  | Human parainfluenza virus 2 | NC_003443.1         | 50         | 38     | 33     |
| Human <del>respirovirus 3</del> | Human parainfluenza virus 3 | EU326526.1          | 50         | 38     | 38     |
| Human <del>rubulavirus 4</del>  | Human parainfluenza virus 4 | NC_021928.1         | 40         | 38     | 33     |
| Influenza A virus               | A/California/07/2009        | NC_026431.1         | 35         | 42     | 29     |
|                                 |                             | NC_026432.1         | 40         | 0      | 29     |
|                                 |                             | NC_026433.1         | 40         | 38     | 29     |
|                                 |                             | NC_026434.1         | 45         | 33     | 0      |
|                                 |                             | NC_026435.1         | 35         | 33     | 33     |
|                                 |                             | NC_026436.1         | 35         | 29     | 38     |
|                                 |                             | NC_026437.1         | 35         | 29     | 33     |
|                                 |                             | NC_026438.1         | 40         | 33     | 33     |

|                                     |                      |               |    |    |    |
|-------------------------------------|----------------------|---------------|----|----|----|
| Influenza B virus                   | B/Lee/1940           | NC_002204.1   | 40 | 38 | 33 |
|                                     |                      | NC_002205.1   | 35 | 50 | 33 |
|                                     |                      | NC_002206.1   | 40 | 29 | 29 |
|                                     |                      | NC_002207.1   | 40 | 33 | 29 |
|                                     |                      | NC_002208.1   | 35 | 33 | 38 |
|                                     |                      | NC_002209.1   | 0  | 33 | 29 |
|                                     |                      | NC_002210.1   | 0  | 33 | 29 |
|                                     |                      | NC_002211.1   | 35 | 29 | 33 |
| Enterovirus D68                     | <del>Fermon</del>    | NC_038308.1   | 40 | 38 | 38 |
| Respiratory syncytial virus         | Strain not indicated | NC_001803.1   | 45 | 38 | 38 |
| Rhinovirus A                        | ATCC VR-1559         | NC_038311.1   | 45 | 38 | 38 |
| <i>Chlamydia pneumoniae</i>         | TW-183               | NC_005043.1   | 60 | 58 | 50 |
| <del>Haemophilus influenzae</del>   | NCTC8143             | NZ_LN831035.1 | 60 | 45 | 50 |
| <i>Legionella pneumophila</i>       | NCTC12273            | NZ_LR134380.1 | 65 | 54 | 67 |
| <i>Mycobacterium tuberculosis</i>   | H37Rv                | NC_018143.2   | 60 | 50 | 50 |
| <del>Streptococcus pneumoniae</del> | NCTC7465             | NZ_LN831051.1 | 65 | 67 | 54 |
| <i>Streptococcus pyogenes</i>       | NCTC8198             | NZ_LN831034.1 | 75 | 54 | 50 |
| <i>Bordetella pertussis</i>         | 18323                | NC_018518.1   | 65 | 50 | 46 |
| <i>Mycoplasma pneumoniae</i>        | FH                   | NZ_CP010546.1 | 60 | 46 | 54 |
| Influenza C virus                   | C/Ann Arbor/1/50     | NC_006312.2   | 0  | 0  | 0  |
|                                     |                      | NC_006311.1   | 0  | 0  | 0  |
|                                     |                      | NC_006310.2   | 0  | 0  | 0  |
|                                     |                      | NC_006309.2   | 45 | 29 | 29 |
|                                     |                      | NC_006308.2   | 45 | 33 | 29 |
|                                     |                      | NC_006307.2   | 40 | 38 | 29 |
|                                     |                      | NC_006306.2   | 40 | 29 | 29 |
|                                     |                      |               |    |    |    |














|                                     |                      |               |    |    |    |
|-------------------------------------|----------------------|---------------|----|----|----|
| <i>Parechovirus</i>                 | Human parechovirus 1 | NC_038319.1   | 40 | 42 | 29 |
| <i>Candida albicans</i>             | SC5314               | NC_032096.1   | 60 | 75 | 50 |
|                                     |                      | NC_032095.1   | 60 | 50 | 50 |
|                                     |                      | NC_032094.1   | 65 | 46 | 50 |
|                                     |                      | NC_032093.1   | 60 | 50 | 50 |
|                                     |                      | NC_032092.1   | 60 | 50 | 50 |
|                                     |                      | NC_032091.1   | 60 | 79 | 50 |
|                                     |                      | NC_032090.1   | 65 | 58 | 63 |
|                                     |                      | NC_032089.1   | 60 | 50 | 54 |
| <i>Corynebacterium diphtheriae</i>  | NCTC11397            | NZ_LN831026.1 | 70 | 54 | 54 |
| <i>Legionella longbeachae</i>       | NSW150               | NC_013861.1   | 70 | 54 | 58 |
| <i>Bacillus anthracis</i> (Anthrax) | Vollum               | NZ_CP007666.1 | 60 | 63 | 50 |
| <i>Moraxella cararrhalis</i>        | BBH18                | NC_014147.1   | 70 | 54 | 63 |
| <i>Neisseria elongata</i>           | ATCC 29315           | NZ_CP007726.1 | 70 | 50 | 50 |
| <i>Neisseria meningitidis</i>       | NCTC10025            | NZ_LR134525.1 | 85 | 54 | 50 |
| <i>Pseudomonas aeruginosa</i>       | PAO1                 | NC_002516.2   | 75 | 54 | 54 |
| <i>Staphylococcus epidermidis</i>   | ATCC 14990           | NZ_CP035288.1 | 60 | 50 | 50 |
| <i>Streptococcus salivarius</i>     | NCTC8618             | NZ_LR134274.1 | 60 | 54 | 54 |
| <i>Leptospira interrogans</i>       | FMAS AW1             | NZ_CP039283.1 | 65 | 54 | 50 |
| <i>Chlamydia psittaci</i>           | 6BC                  | NC_017287.1   | 65 | 46 | 54 |

**Table 10: Microbial Interference Testing.** Ct values are for SARS-CoV-2.

| Organism                      | SARS-CoV-2          | Replicate 1 Ct | Replicate 2 Ct | Replicate 3 Ct | Average Ct |
|-------------------------------|---------------------|----------------|----------------|----------------|------------|
| None                          | 750 GEC/ml (3X LoD) | 33.01          | 33.63          | 33.57          | 33.40      |
| <i>Neisseria meningitidis</i> | 750 GEC/ml (3X LoD) | 33.80          | 33.73          | 34.86          | 34.13      |
| None                          | Negative            | NA*            | NA             | NA             | NA         |
| <i>Neisseria meningitidis</i> | Negative            | NA             | NA             | NA             | NA         |

\* NA: No Amplification

## Symbols Legend

|   |   |   |   |                                  |
|---|---|---|---|----------------------------------|
|  | Collect sample by   |  | See instruction manual  |                                  |
|   |  | Catalog number  |  | Contains sufficient for <n>tests |
|  | Manufacturer  |  | Package contains  |                                  |
|  | Storage instructions  |  | Reagent   |                                  |
|  | Lot number  |  | Positive Control  |                                  |
|  | Do not reuse  |   |   |                                  |
|  | Keep away from sunlight   |   |   |                                  |



IN VITRO DIAGNOSTIC MEDICAL DEVICE

## Technical Support

For technical support, call Zymo Research Corp. Technical Support at 1-949-679-1190 ext. 3, email [tech@zymoresearch.com](mailto:tech@zymoresearch.com)

CFX Maestro is a registered trademark of Bio-Rad Laboratories, Inc.  
KingFisher is a registered trademark of Thermo Fisher Scientific.

Authorized Representative  
Biotik@-Catherine David,  
42 route du périmètre,  
74940 Annecy, France



Zymo Research Corp.  
17062 Murphy Avenue  
Irvine, CA 92614 USA

© 2020 Zymo Research Corp. All Rights Reserved.  
[www.zymoresearch.com](http://www.zymoresearch.com)  
March 2020  
Version 1.0.2

The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®