



ZYMO RESEARCH

*The Beauty of Science is to Make Things Simple®*

# DNase I Set

Recombinant, lyophilized  
RNase-free & protease-free

## Description:

DNase I is an endonuclease that nonspecifically cleaves single- and double-stranded DNA. It requires divalent metal cations to be active.

## Key Features:

- Superior activity.
- Certified RNase-free and protease-free.
- Withstand more than 10 freeze-thaw cycles.

## Storage:

After reconstitution, place on ice until ready to use or store frozen aliquots (-20°C).

## General Applications:

- Remove DNA from protein and RNA samples.
- Preparation of DNA-free RNA.
- DNA labeling by nick-translation.
- Removal of DNA template after in vitro transcription.
- Removal of DNA from RNA samples prior to downstream applications (NextGen sequencing, RT-PCR).

## Components:

Cat.#	Size	Buffer Volume
E1010	250 U	4 ml
E1011	1500 U	4 ml
E1012	1500 U x 5	16 ml

### Unit Definition:

One Kunitz causes an increase in absorbance at 260 nm of 0.001 per minute per ml, at 25°C, pH 7.5, when acting on salmon sperm DNA according to the assay method of Kunitz.

### Reconstitution:

Briefly centrifuge the tube to ensure that the lyophilized DNase I settles at the bottom.

Add DNase/RNase-free water to the lyophilized DNase I, mix by gentle inversion. Avoid phosphate buffer and calcium chelators.

### DNase I Reaction Set Up:

Add DNA digestion buffer (DDB) and DNase I to the sample and incubate at room temperature for 15 min. Make sure the final concentration of DDB is 1x for in-solution reaction. The provided DDB is 10x.

One Unit of DNase I will completely digest 15 µg of genomic DNA in less than 10 minutes at 25°C.

### Inactivation:

Heat inactivate at 75 °C for 10 min with 5 mM EDTA.

Precautions and Disclaimer:

This product is for R&D use only.

### Technical Support:

1-888-882-9682

tech@zymoresearch.com



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple®*

