



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

# Concentration and Nucleic Acid Extraction of Pathogens in Wastewater

Assessment of Compatibility of Ceres Nanosciences' Nanotrap<sup>®</sup> Particles with  
Zymo Research *Quick-DNA/RNA*<sup>™</sup> Water Kit for Wastewater Surveillance

E. Carrasco<sup>1</sup>, L. Saunders<sup>2</sup>, X. Cheng<sup>1</sup>, and K. Locken<sup>1</sup>  
Zymo Research Corporation<sup>1</sup>, Ceres Nanosciences<sup>2</sup>



## INTRODUCTION

The objective of this study is to assess the compatibility of Ceres Nanosciences' Nanotrap® Microbiome Particles and Nanotrap® Enhancement Reagents with Zymo Research's Quick-DNA/RNA™ Water Kit. This integrated workflow is intended to provide researchers and industry professionals with a reliable method for monitoring wastewater, facilitating timely and informed public health interventions.

## METHODS

The study followed the protocol provided by Ceres Nanosciences for using Nanotrap Microbiome A and Microbiome B Particles with Enhancement Reagent 3 (ER3) to concentrate 10 mL of raw wastewater. (Figure 1) The samples were spiked with SARS-CoV-2, Influenza A, *Listeria monocytogenes*, and *Salmonella enterica* to evaluate the recovery efficiency of both viral and bacterial targets.

The concentrated samples were then processed using the Quick-DNA/RNA™ Water Kit, starting with resuspending the Nanotrap Particles with DNA/RNA Shield™ and incubating the samples at 56 °C for 10 minutes. The lysate was subsequently extracted of total DNA/RNA using Zymo-Spin™ technology. The eluted DNA/RNA was subsequently analyzed using digital PCR (dPCR).

## METHOD OVERVIEW

- Concentration Reagents: Nanotrap Microbiome A and B Particles and Enhancement Reagent 3 (ER3)
- Nucleic Acid Extraction Kit: Quick-DNA/RNA™ Water Kit
- Pathogen Spike-In Targets: SARS-CoV-2, Influenza A, *Listeria monocytogenes*, *Salmonella enterica*
- Sample: 10 mL Wastewater from two different sources
- Assay: dPCR in technical duplicate



Figure 1. Nanotrap Microbiome A and B Particles.

## RESULTS

The study demonstrated that Nanotrap Particles are compatible with the Quick-DNA/RNA™ Water Kit. Eluted DNA/RNA was analyzed using dPCR with (assay) on the QIAcuity One – 5plex Digital PCR system.




Figure 2 Viral Detection. Detection of A) SARS-CoV-2 and B) Influenza A in raw wastewater samples by dPCR that were processed using Nanotrap Microbiome Particles and the Quick-DNA/RNA™ Water Kit.

## Viral Recovery

The combined use of Nanotrap Microbiome Particles and the Quick-DNA/RNA™ Water Kit successfully concentrated and extracted viral nucleic acids from wastewater samples. The recovery of SARS-CoV-2 in two wastewater samples was quantified by dPCR at 75.3 and 188.9 copies per microliter (Cp/μL), respectively. (Figure 2) For Influenza A, the recovery was 437.8 and 575.6 Cp/μL, respectively. The variation between the two samples is expected to be due to differences in the two different wastewater samples. These results indicate effective recovery of enveloped, single-stranded RNA viruses CoV-2 representing coronaviruses and Influenza A representing orthomyxoviruses).

## Bacterial Recovery

Similarly, the workflow effectively concentrated and extracted bacterial nucleic acids from the wastewater samples. The recovery of *Listeria monocytogenes* was measured by dPCR at 20.9 and 56.7 Cp/μL, while *Salmonella enterica* was measured at 86.9 and 72.6 Cp/μL. (Figure 3) The variation between the two samples is expected to be due to differences in the two different wastewater samples. These results demonstrate the capability to recover both gram-positive (*Listeria monocytogenes*) and gram-negative bacteria (*Salmonella enterica*), encompassing a range of bacterial types including those that are challenging to lyse.

## CONCLUSION

This study confirms that Ceres Nanosciences' Nanotrap Microbiome Particles are fully compatible with Zymo Research's Quick-DNA/RNA™ Water Kit, establishing an efficient process for concentrating and extracting viral and bacterial nucleic acids from wastewater. The developed workflow offers a powerful tool for wastewater surveillance, enabling accurate detection of pathogens and supporting proactive public health strategies.

## LEARN MORE ABOUT THE QUICK-DNA/RNA™ WATER KIT



[www.zymoresearch.com/pages/quick-dna-rna-water-kit](http://www.zymoresearch.com/pages/quick-dna-rna-water-kit)




Figure 3 Bacterial Detection. PCR detection of A) *L. monocytogenes* and B) *S. enterica* in raw wastewater, utilizing Nanotrap Microbiome Particles and the Quick-DNA/RNA Water Kit for sample processing.