



ZYMO RESEARCH



Quick-RNA™ Tissue/Insect Microprep

RNA from any tough-to-lyse sample

Highlights

- Quick, spin-column purification of total RNA (including small/microRNAs) from any tissue, insect, and arthropod specimens (e.g., mosquitoes, bees, lice, ticks, *Drosophila melanogaster*, etc.)
- ZR BashingBead™ Lysis Tubes are ultra-high density, fracture resistant, chemically inert ceramic beads and used for the robust homogenization of any tough-to-lyse sample.
- RNA is ready for Next-Gen Sequencing, RT/qPCR and any downstream application, etc.

Catalog Numbers:
R2030



Scan with your smart-phone camera to
view the online protocol/video.



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682

Table of Contents

Product Contents	01
Specifications	02
Product Description	03
Protocol	04
(I) Buffer Preparation	04
(II) Total RNA Purification	05
Appendices	06
Samples Stored in DNA/RNA Shield™	06
Protein Purification	07
DNase I Treatment	07
Ordering Information	08
Complete Your Workflow	09
Troubleshooting Guide	10
Notes	11
Guarantee	13

Product Contents

Quick-RNA™ Tissue/Insect Microprep	R2030 (50 prep)
RNA Lysis Buffer	50 ml
RNA Prep Buffer	25 ml
RNA Wash Buffer ¹ (concentrate)	24 ml
DNase/RNase-Free Water	1 ml
ZR BashingBead™ Lysis Tubes (0.1 & 2.0 mm)	50
Zymo-Spin™ IICG Columns	50
Zymo-Spin™ IC Columns	50
Collection Tubes	100
Instruction Manual	1

Storage Temperature - Store all kit components (i.e., buffers, columns) at room temperature.

Before use:

1 Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate.

Specifications

- **Sample Sources** – Up to 10 mg (or $n > 1$) tissue, insect or arthropod specimen (e.g., mosquitoes, bees, lice, ticks, *Drosophila melanogaster*, etc.).
- **Size** – Total RNA including small/microRNAs (≥ 17 nt).
- **Purity** – A_{260}/A_{280} & $A_{260}/A_{230} > 1.8$. RNA is ready for Next-Gen Sequencing, RT/qPCR, etc.
- **Binding Capacity** – **Zymo-Spin™ IC Column** yield up to 10 μ g RNA.
- **Compatibility** – For samples stored in **DNA/RNA Shield™**, see page 6. DNA/RNA Shield™ lyses cells, inactivates nucleases and infectious agents (e.g., virus, pathogens) and is ideal for safe sample storage and transport at ambient temperatures.
- **Elution Volume** – ≥ 6 μ l **DNase/RNase-Free Water**.
- **Equipment Needed (user provided)** – Microcentrifuge, vortex, and a high-speed homogenizer/cell disruptor or bead beater (e.g., MP Bio FastPrep-24, Bertin Precellys, etc.) (recommended).


Product Description

The **Quick-RNA™ Tissue/Insect Microprep Kit** provides for rapid (10-minute) isolation of total RNA from tissue and various insect/arthropod specimens (e.g. mosquitoes, bees, lice, ticks, *Drosophila melanogaster*, etc.). Ultra-high density **ZR BashingBeads™** are included for sample homogenization coupled with a robust buffer system, delivering total RNA (including small/microRNAs).

The **Zymo-Spin™ IICG Column** allows for high-capacity DNA elimination and the subsequent **Zymo-Spin™ IC Column** efficiently binds total RNA.

RNA is eluted in as little as 6 µl and is suitable for subsequent procedures including RT-qPCR.

Efficient Recovery of RNA from any Insect Sample



Total RNA isolation from n=2 *Drosophila sp.* individuals (in duplicate) using the **Quick-RNA™ Tissue/Insect** kit. Samples were processed (2x 30sec at 6 m/s) with FastPrep®-24 Instrument (MP Biomedicals) and resolved along the RNA Millenium™ Markers (Ambion) in a 1% native agarose gel.

Protocol

The protocol consists of: (I) Buffer Preparation and (II) Total RNA Purification

(I) Buffer Preparation

- ✓ Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer**.

(II) Total RNA Purification

- ✓ Perform all steps at room temperature and centrifugation at 10,000-16,000 x g for 30 seconds, unless specified.
- 1. Transfer an appropriate amount (see table) of fresh or frozen tissue into a **ZR BashingBead Lysis Tube** (0.1 & 2.0 mm) and add 800 µl **RNA Lysis Buffer**.

Mammalian Tissue ¹	Insect Tissue	Add RNA Lysis Buffer
High-yield: ≤ 2 mg Low-yield: ≤ 5 mg	≤ 20 mg	800 µl

- 2. Secure tube in a bead beater fitted with a 2 ml tube hold assembly and process. See example below:

Homogenizers	Bead-Beating Time	
	Mammalian	Insect
High-speed (e.g., MP Bio FastPrep-24, Bertin Precellys)	30 – 60 sec	3 – 5 min
Low-speed (e.g., Vortex Genie)	3 – 5 min	15 – 20 min

- 3. Centrifuge the tube for 1 minute to pellet debris.
- 4. Transfer the cleared supernatant into a **Zymo-Spin™ IIICG Column¹** (green) in a **Collection Tube** and centrifuge. Save the flow-through!
- 5. To the flow-through, add an equal volume of ethanol (95-100%) (1:1) and mix well.

Example: Add 400 µl ethanol to 400 µl flow-through.

- 6. Transfer the mixture into a **Zymo-Spin™ IC Column²** in a **Collection Tube** and centrifuge³. Discard the flow-through.

At this point, **DNase I** treatment (optional) can be performed. See page 7.

- 7. Add 400 µl **RNA Prep Buffer** to the column and centrifuge. Discard the flow-through.
- 8. Add 700 µl **RNA Wash Buffer** to the column and centrifuge. Discard the flow-through.
- 9. Add 400 µl **RNA Wash Buffer** and centrifuge the column for 1 minute to ensure complete removal of the wash buffer. Then carefully, transfer the column into a nuclease-free tube (not provided).
- 10. Add 15 µl **DNase/RNase-Free Water** directly to the column matrix and centrifuge. The eluted RNA can be used immediately or stored frozen.

Alternatively, for highly concentrated RNA use ≥ 6 µl elution.

¹ High-yield tissue (e.g., spleen, liver), low-yield tissue (e.g., brain, heart, muscle, lung, intestine, kidney)

² To process samples > 700 µl, columns may be reloaded.

³ At this point, proteins can be purified from the flow-through. See page 7.

Appendices

Samples stabilized and stored in DNA/RNA Shield™

- ✓ If frozen, thaw sample in **DNA/RNA Shield™** to room temperature (20-30°C). Mix well by vortex.

Homogenized Sample

1. Transfer 400 µl of sample homogenized in **DNA/RNA Shield™** to a new RNase-free tube (not provided).
2. Add 400 µl of **RNA Lysis Buffer (1:1)** to the of sample homogenized in **DNA/RNA Shield™** and mix well.
3. Proceed to Total RNA Purification (page 5, step 3).

Non-homogenized Sample

1. Transfer 800 µl - 1 ml of sample suspended in **DNA/RNA Shield™** to a **ZR BashingBead Lysis Tube**.
2. Secure the **ZR BashingBead Lysis Tube** in a bead beater fitted with a 2 ml tube holder assembly and process (see homogenization suggestions on page 5, step 2).
3. Centrifuge the **ZR BashingBead Lysis Tube** for 1 minute at high speed (e.g., 16,000 x g).
4. Transfer 400 µl of the supernatant to a new RNase-free tube (not provided).
5. Add 400 µl of **RNA Lysis Buffer (1:1)** to the supernatant and mix well.
6. Proceed to Total RNA Purification (page 5, step 4).

(Appendices continued)

Protein Purification: Acetone Precipitation of Proteins

- ✓ After the RNA binding to the column (page 5, step 6), the protein content (denatured) in the flow-through can be purified:
- 1. Add 4 volumes of cold acetone (-20°C) to the flow-through (4:1) and mix.
- 2. Incubate the samples for 30 minutes on ice.
- 3. Centrifuge at max speed for 10 minutes. Discard the supernatant. Keep the pellet.
- 4. Add 400 µl ethanol (95-100%) to the protein pellet. Centrifuge at max speed for 1 minute. Discard the supernatant.
- 5. Air-dry the protein pellet for 10 minutes at room temperature.
- 6. Resuspend and vortex the pellet in a buffer appropriate for downstream application (e.g., SDS-PAGE sample loading buffer).

DNase I Treatment (in-column)

- ✓ Perform DNase I treatment with DNase I Set (#E1010) and RNA Wash Buffer (concentrate; #R1003-3-6); available separately.
- 1. Following RNA binding step (page 5, step 6), add 400 µl **RNA Wash Buffer** to the column, centrifuge and discard the flow-through.
- 2. For each sample to be treated, prepare **DNase I Reaction Mix** (see table below) in an RNase-free tube (not provided) and mix by gentle inversion. Then add 40 µl directly into column matrix and incubate at room temperature (20-30°C) for 15 minutes. Proceed with the purification (page 5, step 7).

DNase I Reaction Mix

DNase I (reconstituted; 1 U/µl) ¹	5 µl
DNA Digestion Buffer	35 µl

¹ Prior to use, reconstitute the lyophilized **DNase I** with 275 µl DNase/RNase-Free Water. Mix by gentle inversion and store frozen aliquots. * Unit definition – one unit increases the absorbance of a high molecular weight DNA solution at a rate of 0.001 A₂₆₀ units/ml of reaction mixture at 25°C.

Ordering Information

Product Description	Catalog No.	Size
Quick-RNA™ Tissue/Insect Microprep	R2030	50 preps.

Individual Kit Components	Catalog No.	Amount
RNA Lysis Buffer	R1060-1-50 R1060-1-100	50 ml 100 ml
RNA Prep Buffer	R1060-2-25 R1060-2-100	25 ml 100 ml
RNA Wash Buffer (concentrate)	R1003-3-24 R1003-3-48	24 ml 48 ml
DNase/RNase-Free Water	W1001-10 W1001-30	10 ml 30 ml
DNase I Set (lyophilized) (250 U supplied with DNA Digestion Buffer, 4 ml)	E1010	1 set
Zymo-Spin™ IICG Columns	C1006-50-G	50
Zymo-Spin™ IC Columns	C1004-50	50
Collection Tubes	C1001-50	50
ZR BashingBead™ Lysis Tubes (0.1 & 2.0 mm)	S6014-50	50

Complete Your Workflow

- ✓ For tough-to-lyse samples, use ZR BashingBead Lysis Tubes:

ZR BashingBead Lysis Tubes	
2.0 mm beads #S6003	Plant/animal tissue
0.1 + 0.5 mm beads #S6012	Microbes
0.1 + 2.0 mm beads #S6014	Microbes in tissue/insects

- ✓ For isolation of RNA from any sample:

Quick-RNA kits	
Miniprep Plus #R1057/R1058	$\leq 10^7$ cells, ≤ 50 mg tissue
MagBeads #R2132/R2133	Automatable (Tecan, Hamilton, Kingfisher, etc.)


- ✓ For clean-up (purification) and concentration of any RNA sample. (e.g., from the aqueous phase of TRIzol[®] extractions) or from any enzymatic reaction (e.g., DNase I treated RNA):

RNA Clean & Concentrator kits	
Microprep #R1013-R1014	DNase I Set included
MagBeads #R1082	Automatable (Tecan, Hamilton, Kingfisher, etc.)

- ✓ For NGS:

Zymo-Seq RiboFree Total RNA Library Prep kit	
#R3000	12 preps
#R3003	96 preps

Troubleshooting Guide

Problem	Possible Causes and Suggested Solutions
Precipitation, viscous lysate	<p>Incomplete lysis and/or high-mass input:</p> <ul style="list-style-type: none"> - If precipitation occurs (upon adding ethanol to the lysate) or if the lysate is extremely viscous, increase the volume of DNA/RNA Shield™ and/or RNA Lysis Buffer to ensure complete lysis and homogenization until lysate is transparent (see image). 
Low purity (A_{260}/A_{230} nm, A_{260}/A_{280} nm)	<p>Sample handling:</p> <ul style="list-style-type: none"> - Ethanol and/or salt contamination. After centrifugation steps, carefully remove the column from the collection tube to prevent buffer carryover. Alternatively, blot emptied collection tube with a tissue or towel. - Make sure lysate and wash buffers have passed completely through the matrix of the column. This may require centrifuging at a higher speed and/or longer time. <p>Incomplete lysis and/or cellular debris:</p> <ul style="list-style-type: none"> - Increase the volume DNA/RNA Shield™ and/or RNA Lysis Buffer (proportionally) to ensure complete lysis and homogenization. Be sure to centrifuge any cellular debris and then process the cleared lysate.
Low yield	<p>Sample input:</p> <ul style="list-style-type: none"> - Too much input or incomplete lysis/homogenization can cause cellular debris to clog or overload the column and result in compromised nucleic acid recovery. Use less input material and/or increase the volume DNA/RNA Shield™ and/or RNA Lysis Buffer. <p>High-protein content:</p> <ul style="list-style-type: none"> - Perform Proteinase K treatment to the sample prior to purification. See appropriate sample preparation protocol.
DNA contamination	<p>To remove DNA:</p> <ul style="list-style-type: none"> - Perform in-tube DNase I treatment post-purification, refer to the RNA Clean & Concentrator (Cat. R1013) protocol, page 6, “DNase I treatment before RNA clean-up”. Then, add 150 µl RNA Lysis Buffer to the 50 µl reaction mix (3:1) and mix well. Add an equal volume ethanol (95-100%) (1:1) and mix well. Proceed to purification step 6, page 5. - In the future, Perform in-column DNase I treatment, step 6, page 5. - For future preps, increase the volume of DNA/RNA Shield™ and/or RNA Lysis Buffer to ensure complete lysis and homogenization of the sample.
RNA degradation	<p>To prevent RNA degradation:</p> <ul style="list-style-type: none"> - Immediately collect and lyse fresh sample into DNA/RNA Shield™ and/or RNA Lysis Buffer ensure stability. Homogenized samples can be stored frozen for later processing.

For technical assistance, please contact 1-888-882-9682 or email tech@zymoresearch.com



100% satisfaction guarantee on all Zymo Research products, or your money back.

Zymo Research is committed to simplifying your research with quality products and services. If you are dissatisfied with this product for any reason, please call 1(888) 882-9682.

Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

TM Trademarks of Zymo Research Corporation
Quick-RNA[®] is a registered trademark of Zymo Research Corporation.



ZYMO RESEARCH

The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682