


# Optimize Oral Microbiome Profiling, with Precision


Introducing the new **ZymoBIOMICS® Oral Microbiome Standard**

## Eliminate Bias in Your Workflow

Use the ZymoBIOMICS® Oral Microbiome Standard to assess potential biases in your workflow, with 10 carefully selected oral microbiome-specific microorganisms, designed to challenge every facet of your NGS workflow.



**Figure 1.** The ZymoBIOMICS® Oral Microbiome Standard is accurately quantified to closely match the theoretical profile by shotgun metagenomic analysis.



10 commensal and 2 pathogenic oral bacteria, in staggered abundance




G/C content ranging from 27% - 64%, to challenge amplification bias



4x *Streptococcus* species to challenge species-level identification

## True-to-Life Comparison

Designed using data from real, human oral microbiome samples to emulate a representative oral microbiome profile, as compared to the Human Oral Microbiome Database (HOMD).



## #1 Cited for a Reason

ZymoBIOMICS® standards offer several key advantages over competitors' standards.

Features	ZymoBIOMICS® Oral Microbiome Standard	Supplier A
Abundance Variation	< 15%	< 100%
Foreign Contaminants	< 0.01%	< 10%
Strict Manufacturing QC	✓	✗
BioSafety Level	Non-hazardous	BSL-2 Required
Cells per mL	$1.3 \times 10^9$	$1.2 \times 10^8$



**Scan the QR code to purchase now, or request a FREE sample!**

Product	Cat. No	Size
ZymoBIOMICS® Oral Microbiome Standard	D6332	10 preps



info@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682