


# Purify DNA/RNA in <15 Minutes from any Sample

## Quick-DNA/RNA™ Kit


- **Quick & Easy:** Extract DNA and RNA from the any sample in <15 minutes.
- **Sensitive:** Single cell-level recovery of DNA and RNA.
- **Ultra-Pure:** Ready for Next-generation sequencing, RT-qPCR, arrays, etc.

### Universal Sample Compatibility

The Quick-DNA/RNA™ Plus technology is universal and accommodates any sample input including cultured cells, any tissue, blood, tough-to-lyse samples, etc.




### Highest Yields



DNA and total RNA recovery is higher using the Quick-DNA/RNA™ Microprep Plus Kit compared to a Supplier Q kit. Nucleic acids were extracted from 50K HeLa cells (n=2).

### Single-Cell Detection



HeLa cells were titrated down to a single cell, and DNA and total RNA were extracted using the Quick-DNA/RNA™ Microprep Plus Kit. Analysis by RT-qPCR shows high linear recovery of DNA & RNA down to the single-cell level (n=2).

Product	Cat. No.	Size	Binding Capacity	Minimum Elution	Sample Input
Quick-DNA/RNA™ Microprep Plus Kit	D7005*	50 preps	10 µg	6 µl	Cells, Soft Tissue
Quick-DNA/RNA™ Miniprep Kit	D7001	50 preps	25 µg	25 µl	
Quick-DNA/RNA™ Miniprep Plus Kit	D7003T D7003	10 preps 50 preps	100 µg	50 µl	Cells, Any Tissue, Whole Blood
Quick-DNA/RNA™ Magbead Kit	R2130 R2131	1 x 96 preps 4 x 96 preps	20 µg	50 µl	

