


Monkeypox virus (MPXV) inactivation with Zymo Research's DNA/RNA Shield™



nd = not detectable; LOD = limit of detection due to toxicity of buffers

Viral plaque assay demonstrating that DNA/RNA shield effectively inactivates Monkeypox Virus (below the assay limit of detection). Zymo Research, in collaboration with the Institute of Virology, University Medical Center Freiburg, Germany.

Product	Cat.	No. Size
DNA/RNA Shield™	R1100-250	250 ml
DNA/RNA Shield™ (2X Concentrate)	R1200-125	125 ml
DNA/RNA Shield™ SafeCollect™ Swab Collection Kit, 2ml	R1161	1 collection kit