



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

PRODUCT INFORMATION

Isopropyl- β -D-thiogalactopyranoside (IPTG) *

Catalog Nos. **I1001-5 & I1001-25**

Molecular Formula: C₉H₁₈O₅S

Molecular Weight: 238.3

CAS Number: 367-93-1

Concentration: 0.5M in water

Storage: -20 °C

Purity: $\geq$ 98%

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

