



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple™

ZR small-RNA Ladder

Catalog No:

R1090

Content:

ZR small-RNA Ladder, 10 µg

Description:


The ZR small-RNA Ladder is a microRNA size marker for use in polyacrylamide gel separation methods and small RNA size approximation.

Specifications:

The ladder consists of four single-stranded RNA oligo nucleotides: 17, 21, 25 and 29 bases in length.

Storage Temp.:

Upon arrival, store aliquots at -20°C to avoid freeze-thaw cycles.



350 ng **ZR small-RNA Ladder** resolved on a 25% native PAGE gel and stained with GelStar® for 5 minutes.