



ZR DNA Sequencing Clean-up Kits™

- **Long Read Lengths:** Complete elimination of “dye blobs” for high-quality Phred scores.
- **Reusable:** Columns can be reused.
- **Flexible Elutions:** 6-20 µl elution volumes allow for direct loading of samples without precipitation or drying.

Simple, Innovative Workflow



Ultra-pure DNA without “Dye Blobs”



Sequencing chromatogram of pGEM® DNA generated using an ABI 3730xl DNA analyzer. DNA was labeled with ABI BigDye® v3.1 Terminators and cleaned using the ZR DNA Sequencing Clean-up Kit™.

Product	Cat. No.	Size
ZR DNA Sequencing Clean-Up Kits™	D4050	50 preps
	D4051	200 preps
ZR-96 DNA Sequencing Clean-Up Kits™	D4052	2 x 96 preps
	D4053	4 x 96 preps