


GET IT DIRECT

DNA & RNA from TRIzol® in 7 minutes

Direct-zol™ DNA/RNA (any sample: blood, cells, tissue, etc.)

- **One Input, One Column:** Purify DNA & RNA directly from TRIzol® with just one spin-column.
- **Easy Handling:** No phase-separation or precipitation steps.
- **High-Quality:** DNA & RNA (including small & micro RNA) ready for Next-Gen Sequencing, qPCR and RT-qPCR, hybridization, etc.

NGS-Ready DNA and RNA



High quality DNA and RNA purified in duplicate from the same input of mammalian (HeLa) and bacterial (*E. coli*) cells using the **Direct-zol™ DNA/RNA Miniprep Kit**. Samples were visualized using the Agilent 2200 TapeStation® system.

DNA and RNA *directly* from TRIzol®



in separate fractions

Product	Cat. No.	Size	Binding Capacity	Minimum Elution
Direct-zol™ DNA/RNA Miniprep Kit	R2080T R2080 R2081*	10 preps 50 preps 50 preps	25 µg DNA and 50 µg RNA	25 µl

*Supplied with TRI Reagent®