

DNA/RNA Shield™

SafeCollect Saliva Collection Kit

Cat. No. R1211

- DNA/RNA Shield™ will stabilize the RNA content of samples for up to 30 days at ambient temperatures.
- Samples stored in DNA/RNA Shield™ can be frozen at -20°C or -80°C for long-term storage.
- Samples stored in DNA/RNA Shield™ are compatible with most extraction workflows.

Recommended Purification Kits

	Host Cells	Viral	Microbial
DNA	Quick-DNA™ Plus Kits D4068 (Column) D4081 (Magbead)	Quick-DNA™ Viral Kit D3015	ZymoBIOMICS™ DNA Kits D4300 (Column) D4302 (Magbead)
RNA	Quick-RNA™ Kits R1057 (Column) ¹ R2132 (Magbead) ²	Quick-RNA™ Viral Kit R1034	ZymoBIOMICS™ RNA Kits R2001 (Column) R2137 (Magbead)
DNA/RNA	Quick-DNA/RNA™ Kits D7003 (Column) R2130 (Magbead)	Quick-DNA/RNA™ Viral Kit D7020 (Column) R2140 (Magbead)	ZymoBIOMICS™ DNA/ RNA Miniprep Kits R2002 (Column) R2135 (Magbead)

1. Following the "Whole-Blood" extraction method is recommended. The Quick-RNA™ Whole Blood kit can be used as well.
2. It is recommended to centrifuge the sample prior to extraction to remove unlysed bacteria from the sample.

For technical inquiries please contact tech@zymoresearch.com

