




ZYMO RESEARCH

The Beauty of Science is to Make Things Simple



Laboratory Automation Made Simple

Innovation. Pure & Simple.

Our Promise to You

Since its inception in 1994, Zymo Research has been proudly serving the scientific community by providing innovative, reliable, and high quality research tools at affordable prices. Our vision “*The **Beauty of Science** is to Make Things **Simple***” is now truer than ever. Whether it’s epigenetics, microbiomics, DNA, RNA, *E. coli*, or yeast based research, our philosophy remains the same: To provide the highest quality products in the industry while ensuring they are both simple to use and reliable in their performance.

Zymo Research stands on three pillars which form the foundation of our company: Innovation, Quality, and Customer Service. These pillars are fundamental to our culture and ensure our products meet your needs.

Innovation

Zymo Research is globally recognized for its innovation of high quality nucleic acid purification technologies. Some key innovations include streamlining DNA methylation detection, pioneering the micro-elution column for DNA and RNA purification, developing the simplest and the most sophisticated methods for high-quality plasmid DNA purification, and patenting the first RNA purification directly from TRIzol® without phase separation.

Quality

We are committed to quality and guarantee that all of our products and services will meet and exceed your expectations.



Customer Service

We strive for excellence in supporting your scientific endeavors. We pledge to be honest and responsible and treat you as we would like to be treated. Together, we will build a brighter future.

Sustainable Solutions

“Zymo Research is engineering sustainable solutions in the biomedical and biotechnology world. We are committed to making fundamental changes to eliminate cold-chain logistics for biological sample transportations globally, one of the most environmentally contaminating practices in the industry. Through scientific innovation, we aim to reduce environmental waste and create a more sustainable future.” – Dr. Jia, CEO & Founder of Zymo Research




Contents



Automation Support at Zymo Research

1 - 4



Featured Automation Workflows

Fecal Microbial Analysis for Precision Medicine
Early Disease Detection from Cell-Free DNA
Epigenetic Biomarker Discovery in Breast Cancer Research

5 - 10



Custom Manufactured Automation Solutions

11




Zymo X OpenTrons - The First Multiomic Workstation

12



Load N' Go - Pre-filled Automation Kits

13



All Automation Products

Microbiome
Epigenetics & Methylation
Plasmid
Liquid Biopsy & Cell-Free
RNA from TRIZOL
General Nucleic Acid Purification

14 - 17

We Can Help with Laboratory Automation

Automating a laboratory workflow process comes with many challenges:

- Steep Learning Curve
- Time Investment
- Trial & Error

Zymo Research eliminates these hurdles with ready-to-use solutions so you can focus on the science, not the setup.



Free Consultations - Find Your Complete Solution

Collaborate with the experts at Zymo Research to find your automated solution.



Free Scripts on All Popular Open Platforms

Take advantage of our pre-made scripts available for Opentrons FLEX™, Accuris IsoPure™ 96, Hamilton STAR™, KingFisher™ Flex, Tecan Fluent®, and More!



Free Kits for Verification - Try Before You Buy

Want to try before you buy? Validate your workflow at no cost to guarantee performance before commitment.

“My questions are always answered promptly and it is much appreciated. The people providing technical support are always very helpful and friendly. Thank you very much!”

- Jason V. at Mayo Clinic



Consultations for Complete Solutions



Explore Verified Processes



Personalized Support



Complete End-to-End Workflow

1

Select Your Sample Type, Analyte, Application, & Instrument

Scan to Get Started!



The form includes dropdown menus for 'Sample Type' (with options: Animal, Environmental, Human, Microbial), 'Analyte', 'Application', and 'Instrument'. It also features a 'Reset' button and a yellow 'Show Workflow' button. A hand cursor icon is shown over the 'Sample Type' dropdown.

2

Find the Products You Need from Sample Collection to Downstream Application

The product catalog is organized into three sections: 'Collection', 'Extraction/Purification', and 'Application Kit'.
- **Collection:** DNA/RNA Shield Fecal Collection Tube (Cat#: R1101), 10 pack.
- **Extraction/Purification:** ZymoBIOMICS 96 MagBead DNA Kit (Lysis Rack) (Cat#: D4302), 2 x 96 Preps (Lysis Rack) and 2 x 96 Preps (Lysis Tubes).
- **Application Kit:** Quick-16S Plus NGS Library Prep Kit (V1-V2) with Primer Set 1 (Cat#: D6434-PS1), 96 rxns (V1-V2). Below this are four smaller product cards for different versions of the Quick-16S Plus NGS Library Prep Kit (V1-V2, V1-V3, V3-V4, UDI, and V4).
A yellow 'Request a Quote' button is located at the bottom.

3

This Generates Your Automation Workflow!



Schedule a Personalized Consultation Today!

Or Contact: Automation@ZymoResearch.com

Free Scripts for Your Instrument



Save Time on Programming



Zero-Cost Scripting



All Open Platforms Supported

Save time and resources with our free, ready-made scripts for popular automation platforms like Hamilton, Tecan, KingFisher, and more! Seamlessly integrate Zymo assays into your workflow without the need for custom script development.

Have unique parameters for your workflow? Reach out to one of our specialists and we can work to adapt the script to your needs at no cost.

“The support and information I received was great and helped us reduce time spent in feasibility testing and development.”

- Kevin K. at Exact Sciences



Get Your Scripts Today!

Or Contact: Automation@ZymoResearch.com



Free Kits for Trial

Validate your automation workflow with free sample kits, offering a risk-free way to test and optimize your workflow.

Free access to sample kits allows you to easily evaluate and refine your workflow, minimizing upfront investment, streamlining implementation, and ensuring a smooth transition to automation in your lab.



**No-Cost
Verification**



**Complete
Support**



**Optimize
Before You Buy**

“They answered all my automation questions in a prompt manner. They are an excellent resource which helped me implement automation using Zymo kits. Thanks for your help and support!”

- Paulo L. at Stanford Medicine



Fecal Metatranscriptomic Analysis for Precision Medicine



Advance precision medicine with reliable microbiome metatranscriptomic profiling to uncover microbial biomarkers. Accelerate clinical trials and therapeutic discoveries with our end-to-end automated solution, featuring DNA/RNA Shield™ for sample stabilization, ZymoBIOMICS MagBead Kits for complete microbial representation, and Zymo-Seq RiboFree Total RNA Library Kit for RNA-Seq.



1 Fecal Collection

Collecting fecal samples is simple and reliable with the **Bunny Wipe with DNA/RNA Shield™ Fecal Collection Tube (R1138)**.

The Bunny Wipe provides an easy-to-use sampling method, while DNA/RNA Shield™ Reagent preserves the integrity of microbial nucleic acids at room temperature, ensuring high-quality samples for downstream metatranscriptomic analysis.

2 RNA Purification


The **ZymoBIOMICS MagBead RNA Kit (R2137/R2138)** enables efficient, bias-free RNA extraction from fecal samples, while the **RNA Clean & Concentrator MagBead Kit (R1081/R1082)** adds an additional layer of insurance for purified, high-quality RNA, ready for accurate metatranscriptomic analysis. Together, they provide a streamlined, automated extraction workflow that ensures quality input for library prep.

3 Library Prep for NGS

The **Zymo-Seq RiboFree Total RNA Library Kit (R3000)** delivers high-quality, ribosomal RNA-depleted libraries, enabling comprehensive metatranscriptomic analysis of microbial activity and functional pathways while being fully automatable for increased efficiency.


**ZymoBIOMICS RNA Has Over
3,500 Citations!**

Automated on the Opentrons Flex™




Metatranscriptomic Analysis for Precision Medicine Applications.

Functional Pathway Analysis Heatmap with Clustering



This heatmap showcases the abundance of key microbial pathways (e.g., methanogenesis, glycolysis, amino acid biosynthesis) across multiple fecal samples. With this automated metatranscriptomics workflow from Zymo Research, active metabolic functions can be uncovered at unprecedented speed and scale for personalized healthcare and wellness.

Early Disease Detection from Cell-Free DNA



Streamline cfDNA analysis for liquid biopsy applications. Using plasma derived from industry-standard Streck Cell-Free DNA BCT® tubes, the MAGicBead™ cfDNA Kit enables efficient automated cfDNA extraction, while the Zymo-Seq Trio WGBS Kit prepares high-quality libraries for Bisulfite NGS, facilitating precise biomarker detection for early disease detection.

1 Plasma Collection

Blood plasma derived from **Streck Cell-Free DNA BCT® tubes**, the industry standard for liquid biopsy sample collection, is seamlessly compatible with Zymo's MAGicBead cfDNA Kit, ensuring efficient and reliable cfDNA extraction for downstream analysis.

2 cfDNA Extraction

The **MAGicBead™ cfDNA Kit (D4086)** offers superior yield and the fastest extraction protocol in the industry. Its optimized, proprietary magnetic bead-based system ensures maximum recovery of high-quality cfDNA, enabling efficient and reliable preparation for downstream applications.

3 Library Prep for WGBS

The **Zymo-Seq Trio WGBS Kit (D5462/D5463)** enables comprehensive whole-genome bisulfite sequencing (WGBS) of cell-free DNA, providing high-resolution methylation profiles for liquid biopsy applications. By accurately detecting both hypermethylated and hypomethylated regions, the kit facilitates the discovery and analysis of epigenetic biomarkers for early disease detection, monitoring, and classification.




Zymo cfDNA Extraction has over 2,400 Citations!

Automated on the KingFisher™ Flex




Differential Methylation Profiles from cfDNA for Liquid Biopsy



This heatmap illustrates the differential DNA methylation patterns between healthy and cancer samples, grouped by methylation status (hypermethylated or hypomethylated regions). Each row represents a genomic region, while the columns correspond to individual samples. The clustering patterns highlight distinct methylation profiles, with cancer samples showing notable deviations compared to healthy controls, emphasizing the potential of methylation biomarkers for disease detection and classification.

Biomarker Discovery in Breast Cancer Research



Unlock epigenetic insights from challenging FFPE breast tissue samples with ease. Our end-to-end automated workflow featuring our EZ DNA Methylation Lightning Automation Kit delivers reproducible and reliable methylation profiles that are NGS-ready for the Illumina® Infinium™ MethylationEPIC BeadChip Kit.



1

Sample Processing

Zymo Research's **Deparaffinization Solution (D3067-1-20)** and **2X Digestion Buffer (D3050-1)** simplify FFPE tissue processing by efficiently removing paraffin and preparing samples for nucleic acid extraction, ensuring high-quality input for epigenetic workflows.

2

Nucleic Acid Extraction

The **Quick-DNA/RNA MagBead Kit (R2130/R2131)** enables efficient, high-throughput nucleic acid extraction from FFPE samples, providing high-quality fractional DNA and RNA from the same sample, allowing for immediate methylation analysis without additional purification.

3

DNA Methylation

The **EZ DNA Methylation-Lightning Automation Kit (D5049)** streamlines bisulfite conversion, delivering fast, reliable, and high-throughput preparation of DNA for precise methylation analysis in automated epigenetic workflows.


**EZ DNA Methylation-Lightning Has Over
1,400 Citations!**

Automated on the Hamilton Microlab® STAR™



Ready for Profiling with the Illumina® Infinium™ MethylationEPIC BeadChip Kit

Reproducible Automation Results



With an R^2 value of 0.992, this data highlights the EZ DNA Methylation-Lightning Automation Kit's precision and reliability for consistent results, even in high-throughput workflows. Perfect for studies demanding accuracy, such as biomarker discovery in breast cancer research.

Tailored Automation Solutions



Custom Automated Workflows

Not all science fits in a kit. Custom scripts & protocols are available for your unique needs.



Bulk Reagent Volumes

Zymo Research offers **ISO 9001:2015** manufacturing for custom, large scale solutions.



Collaborative Support

Direct access to in-house scientists & automation experts to help you achieve your goals.



Custom Pre-Filled Formats Available



Custom OEM Labels and Packaging Available



Compatible with All Open Deck Platforms



Available in 96 and 16 prep sizes



Schedule to Start Today!

Or Contact: Automation@ZymoResearch.com

Experience Freedom with Opentrons: The Multiomic Workstation.

Zymo Research has partnered with Opentrons to deliver a turn-key, multiomic workstation that provides NGS libraries ready for both metagenomics and metatranscriptomics - from the same sample.



Isolate Fractional DNA & RNA from the same source.

A truly **multiomic solution** that enables metagenomic and metatranscriptomic applications from the same sample.



A Complete Workflow from Sample to Library.

An **end-to-end workflow** that handles sample processing, library prep, and everything in between.



A Turn-Key, Standardized Solution.




Pre-filled **Load N' Go** reagent plates and **pre-scripted** workflows enable a truly standardized solution that works out of the box.



Find Out if it's Right for You!

Or Contact: Automation@ZymoResearch.com

Simplify & Standardize with the Load N' Go Format

-  Standardized, Turn-key Extraction
-  Free Scripts Available
-  Save Time & Money




The Load N' Go Advantage

The Pre-filled Load N' Go format eliminates repetitive labor and error potential for high-throughput automated workflows while standardizing the process out of the box.

Consolidate consumable usage to just 1 rack of tips for both liquid handling and magnetic bead transfer platforms.

Hands on Time Comparison



Note: Figures are from Plate Preparation of the Quick-DNA/RNA HT and the Load N' Go Quick-DNA/RNA HT kits for the KingFisher™ Flex.

Compatible with All Open Platforms

Compatible with Opentrons FLEX™, Accuris IsoPure™ 96, Hamilton STAR™, KingFisher™ Flex, Tecan Fluent®, and More!



| Field of Interest | Kit Name | Cat. Number | Analyte | Prep Size |
|---------------------------|--|-------------|---------|-----------|
| Nucleic Acid Purification | Load N' Go Quick-DNA/RNA HT | R2152 | DNA/RNA | 96 |
| Microbiome | Load N' Go ZymoBIOMICS DNA MagBead Kit | D4311 | DNA | 96 |
| Viral | Load N' Go Quick-DNA/RNA Viral Kit | R2143 | DNA/RNA | 96 |

All Automation Products

Microbiome



| Sample Collection & Preservation | Cat. No. | Size |
|---|-----------------------|-------------------------|
| DNA/RNA Shield™ | R1100-50 | 50 mL |
| DNA/RNA Shield™ (2x Concentrate) | R1200-25 | 25 mL |
| DNA/RNA Shield™ Fecal Collection Tube | R1101 | 10 Pack |
| Bunny Wipe Fecal Sample Collector | R1133-10 | 10 Wipes |
| DNA/RNA Shield™ SafeCollect Swab Collection Kit | R1161 | 1 Collection Kit |
| DNA/RNA Shield™ Collection Tube w/ Swab | R1107 | 50 Pack |
| DNA/RNA Shield™ SafeCollect Saliva Collection Kit | R1211 | 1 Collection Kit |
| DNA/RNA Shield™ Blood Collection Tube | R1150 | 50 Pack |
| Sample Processing & Lysis | Cat. No. | Size |
| ZR BashingBead Lysis Tubes (0.1 & 0.5 mm) | S6012-50 | 50 Pack |
| ZymoBIOMICS BashingBead Lysis Rack (0.1 & 0.5 mm) | S6002-96-6 | 1 Rack |
| ZymoBIOMICS BashingBead Lysis Rack (0.1 & 0.5 mm) w/ Heat Sealing Foils | S6002-96-7 | 1 Rack |
| Extraction Kits | Cat. No. | Size |
| ZymoBIOMICS 96 MagBead DNA Kit | D4302 | 2 x 96 Preps |
| ZymoBIOMICS 96 MagBead DNA Kit - Dx | D4302-E | 2 x 96 Preps |
| Load N' Go ZymoBIOMICS DNA MagBead Kit | D4311 | 96 Preps |
| ZymoBIOMICS MagBead RNA | R2137/R2138 | 96 Preps / 4 x 96 Preps |
| ZymoBIOMICS MagBead DNA/RNA | R2135/R2136 | 96 Preps / 4 x 96 Preps |
| Library Prep | Cat. No. | Size |
| Quick-16S Plus NGS Library Prep Kit (v1-v3) | D6440-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v1-v2) | D6434-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v3-v4, UDI) | D6421-PS1 - D6421-PS8 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v4) | D6430-PS1 - D6430-PS4 | 96 Preps |
| Quick-ITS Plus NGS Library Prep Kit (UDI) | D6425-PS1 - D6425-PS4 | 96 Preps |
| Controls | Cat. No. | Size |
| ZymoBIOMICS Microbial Community Standard | D6300 | 10 Preps |
| ZymoBIOMICS Fecal Reference with TruMatrix™ Technology | D6323 | 10 Preps |
| ZymoBIOMICS Gut Microbiome Standard | D6331 | 10 Preps |
| ZymoBIOMICS Oral Microbiome Standard | D6332 | 10 Preps |

Epigenetics & DNA Methylation



| Sample Collection & Preservation | Cat. No. | Size |
|---|----------|------------------|
| DNA/RNA Shield™ | R1100-50 | 50 mL |
| DNA/RNA Shield™ (2x Concentrate) | R1100-25 | 25 mL |
| DNA/RNA Shield™ SafeCollect Swab Collection Kit | R1161 | 1 Collection Kit |
| DNA/RNA Shield™ Collection Tube w/ Swab | R1107 | 50 Pack |
| DNA/RNA Shield™ SafeCollect Saliva Collection Kit | R1211 | 1 Collection Kit |
| DNA/RNA Shield™ Blood Collection Tube | R1150 | 50 Pack |

| Sample Processing & Lysis | Cat. No. | Size |
|---|------------|---------|
| ZR BashingBead Lysis Tubes (0.1 & 0.5 mm) | S6012-50 | 50 Pack |
| ZR BashingBead Lysis Tubes (0.1 & 2.0 mm) | S6014-50 | 50 Pack |
| Deparaffinization Solution | D3067-1-20 | 20 mL |
| 2X Digestion Buffer | D3001-2-20 | 20 mL |
| Proteinase K w/ Storage Buffer Set | D3001-2-60 | 60 mg |

| Extraction Kits | Cat. No. | Size |
|----------------------------|-------------|-------------------------|
| Quick-DNA Magbead Plus Kit | D4081/D4082 | 96 Preps / 4 x 96 Preps |
| Quick-DNA/RNA Magbead Kit | R2130/R2131 | 96 Preps / 4 x 96 Preps |

| Methylation Kits | Cat. No. | Size |
|---|-------------|-----------------------------|
| EZ-96 DNA Methylation-Lightning MagPrep | D5046/D5047 | 4 x 96 Preps / 8 x 96 Preps |
| EZ-96 DNA Methylation MagPrep | D5040/D5041 | 4 x 96 Preps / 8 x 96 Preps |
| EZ-96 DNA Methylation-Gold MagPrep | D5042/D5043 | 4 x 96 Preps / 8 x 96 Preps |
| EZ-96 DNA Methylation-Direct MagPrep | D5044/D5045 | 4 x 96 Preps / 8 x 96 Preps |

| Library Prep | Cat. No. | Size |
|----------------------------------|-------------|---------------------|
| Zymo-Seq Trio WGBS Library Kit | D5462/D5463 | 24 Preps / 96 Preps |
| Zymo-Seq WGBS Library Kit | D5465 | 24 Preps |
| Zymo-Seq RRBS Library Prep Kit | D5460/D5461 | 24 Preps / 48 Preps |
| Zymo-Seq ATAC Library Kit | D5458 | 12 Preps |
| Pico Methyl-Seq Library Prep Kit | D5455/D5456 | 10 Preps / 25 Preps |
| RRHP 5-hmC Library Prep Kit | D5450/D5451 | 12 Preps / 25 Preps |

Plasmid



| Extraction Kits | Cat. No. | Size |
|---------------------------------------|-------------|-------------------------|
| Zyppy-96 Plasmid MagBead Miniprep Kit | D4100/D4101 | 96 Preps / 4 x 96 Preps |

Liquid Biopsy & Cell-Free



| Extraction Kits | Cat. No. | Size |
|--------------------------------|----------|-----------------|
| MAGicBead™ cfDNA Isolation Kit | D4086 | 2 mL x 50 Preps |

| Library Prep | Cat. No. | Size |
|--------------------------------|-------------|---------------------|
| Zymo-Seq Trio WGBS Library Kit | D5462/D5463 | 24 Preps / 96 Preps |

RNA from TRIZOL



| Sample Collection & Preservation | Cat. No. | Size |
|----------------------------------|-------------|--------|
| TRI Reagent | R2050-1-200 | 200 mL |

| Extraction Kits | Cat. No. | Size |
|---------------------------|-------------|-------------------------|
| Direct-zol-96 MagBead RNA | R2100/R2102 | 96 Preps / 4 x 96 Preps |

| Library Prep & PCR | Cat. No. | Size |
|---|-------------|---------------------|
| Zymo-Seq RiboFree Total RNA Library Kit | R3000/R3003 | 12 Preps / 96 Preps |
| ZymoScript One-Step RT-qPCR Kit | R3014 | 100 Reactions |

General Nucleic Acid Purification



| Sample Collection & Preservation | Cat. No. | Size |
|---|----------|------------------|
| DNA/RNA Shield™ | R1100-50 | 50 mL |
| DNA/RNA Shield™ (2x Concentrate) | R1200-25 | 25 mL |
| DNA/RNA Shield™ Blood Collection Tube | R1150 | 50 Pack |
| DNA/RNA Shield™ SafeCollect Swab Collection Kit | R1161 | 1 Collection Kit |
| DNA/RNA Shield™ Collection Tube w/ Swab | R1107 | 50 Pack |
| DNA/RNA Shield™ SafeCollect Saliva Collection Kit | R1211 | 1 Collection Kit |

| Extraction Kits | Cat. No. | Size |
|-----------------------------|-------------|-------------------------|
| Quick-DNA Magbead Plus | D4081/D4082 | 96 Preps / 4 x 96 Preps |
| Quick-DNA HMW Magbead | D6060 | 96 Preps |
| Quick-RNA MagBead | R2132/R2133 | 96 Preps / 4 x 96 Preps |
| Quick-DNA/RNA MagBead | R2130/R2131 | 96 Preps / 4 x 96 Preps |
| Quick-DNA/RNA HT | R2150/R2151 | 250 Preps / 1,000 Preps |
| Quick-DNA/RNA Viral MagBead | R2140/R2141 | 250 Preps / 1,000 Preps |

| Library Prep & PCR | Cat. No. | Size |
|--|-----------------------|----------------------|
| Zymo-Seq Trio WGBS Library Kit | D5462/D5463 | 24 Preps / 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v1-v3) | D6440-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v1-v2) | D6434-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v3-v4, UDI) | D6421-PS1 - D6421-PS8 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v4) | D6430-PS1 - D6430-PS4 | 96 Preps |
| Quick-ITS Plus NGS Library Prep Kit | D6424-PS1 - D6424-PS4 | 96 Preps |
| Zymo-Seq RiboFree Total RNA Library Kit | R3000/R3003 | 12 Preps / 96 Preps |
| ZymoTaq DNA Polymerase | E2001/E2002 | 50 Preps / 200 Preps |
| ZymoScript One-Step RT-qPCR Kit | R3014 | 100 Reactions |

Clean-up & Library Prep



| Clean-Up / Purification | Cat. No. | Size |
|--|----------|-----------|
| DNA Clean & Concentrator MagBead Kit | D4012 | 100 Preps |
| RNA Clean & Concentrator MagBead Kit | R1081 | 96 Preps |
| RNA Clean & Concentrator MagBead Kit (Supplied with DNase I Set) | R1082 | 96 Preps |

| Library Prep & PCR | Cat. No. | Size |
|--|-----------------------|---------------------|
| Zymo-Seq RiboFree Total RNA Library Kit | R3000/R3003 | 12 Preps / 96 Preps |
| ZymoScript One-Step RT-qPCR Kit | R3014 | 100 Reactions |
| Zymo-Seq WGBS Library Kit | D5465 | 24 Preps |
| Zymo-Seq RRBS Library Prep Kit | D5460/D5461 | 24 Preps / 48 Preps |
| Zymo-Seq ATAC Library Kit | D5458 | 12 Preps |
| Pico Methyl-Seq Library Prep Kit | D5455/D5456 | 10 Preps / 25 Preps |
| RRHP 5-hmC Library Prep Kit | D5450/D5451 | 12 Preps / 25 Preps |
| Zymo-Seq Trio WGBS Library Kit | D5462/D5463 | 24 Preps / 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v1-v3) | D6440-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v1-v2) | D6434-PS1 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v3-v4, UDI) | D6421-PS1 - D6421-PS8 | 96 Preps |
| Quick-16S Plus NGS Library Prep Kit (v4) | D6430-PS1 - D6430-PS4 | 96 Preps |
| Quick-ITS Plus NGS Library Prep Kit | D6424-PS1 - D6424-PS4 | 96 Preps |



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple®

*The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®*



info@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682