

DNA/RNA Co-Purification

Product Chart

Specifications	Quick-DNA/RNA™ Miniprep Kit & Microprep Kit	sDNA/RNA Clean & Concentrator™	Quick-DNA/RNA™ Viral Kit	Quick-DNA/RNA™ Viral 96 Kit
Format	Spin-Column	Spin-Column	Spin-Column	96-Well
Binding Capacity	25 µg DNA 25 µg RNA	10 µg	25 µg	10 µg
Elution Volume	≥ 50 µl DNA ≥ 25 µl RNA	≥ 6 µl	≥ 35 µl	≥ 10 µl
Processing Time	15 min.	10 min.	5 min.	15 min.
Applications				
Parallel Purification	•			
Co-Purification		•	•	•
Fresh/Frozen Soft Tissue	•			
Fresh/Frozen Solid Tissue	limited			
Bacteria	limited			
Yeast	limited			
Buffy Coat	•			
Cultured Cells	•			
Small RNA	•	•		
Probe Purification		•		
Whole Blood (≤ 50 µl)			•	•
Plasma/Serum			•	•
Virus			•	•
Page Number	95	96	97	97

New Products:
 Quick-DNA/RNA Viral MagBead Kit
 Quick-DNA/RNA Pathogen MagBead Kit
 Quick-cfDNA/RNA Serum Plasma Kit