

DNA/RNA Shield Fecal Collection Tube

Cat no(s): R1101

All biological samples stabilized in DNA/RNA Shield™ are:

- (a) compatible with most workflows (e.g., Qiagen, Roche, Thermo-Fisher).
- (b) a direct input for high-throughput processing (e.g., Hamilton, Tecan, bioMérieux, PerkinElmer, Eppendorf, Promega, Opentrons).

Recommended Purification kits

DNA	ZymoBIOMICS DNA Kits D4300 (Column) D4302 (Magbead)
RNA	ZymoBIOMICS RNA Miniprep Kit R2001
DNA/RNA	ZymoBIOMICS DNA/RNA Miniprep Kit R2002

For technical inquiries please contact tech@zymoresearch.com

