



**ZYMO RESEARCH**

*The Beauty of Science is to Make Things Simple*

## Product Change Notification Letter

**Date:** February 21, 2025

**Product Change Notification #:** RD.PCN.134

Dear Valued Distributor/Customer:

This letter serves to inform you that several of our environmental DNA and RNA extraction kits will undergo a change of the Zymo-Spin III-HRC Filters (C1058-50) for the new updated version Zymo-Spin IV-IR HRC Filters (C1010-50). As part of our commitment to quality and continuous process improvements, we have integrated customer feedback and recommendations into the new design. The revised inhibitor removal columns feature a faster workflow, improved appearance, and more robust inhibitor removal performance.

What will change with this update:



- ✓ **Simplified Workflow and Improved Appearance:** HRC columns will now come with a slurry or “wet” format, deleting the need to add prep buffer prior to adding your sample.
- ✓ **Improved Performance:** Updated chemical formulation allows the HRC column to bind to more PCR inhibitors effectively cleaning up strongly inhibited samples in just one step!

We will be transitioning to the new updated version beginning with the **Quick–DNA Plant/Seed Miniprep kit (Cat # D6020)** on **June 1<sup>st</sup>, 2025**. The Zymo-Spin III-HRC Filters (C1058-50) will still be available for purchase until July of 2025. If you need to complete your study using the Quick–DNA Plant/Seed Miniprep kit (Cat # D6020) using the Zymo-Spin III-HRC Filter format, please purchase prior to June 1<sup>st</sup>, 2025.

17062 Murphy Ave. • Irvine, CA 92614 • USA

Tel: 1-888-882-9682 • Fax: 1-949-266-9452 • [www.zymoresearch.com](http://www.zymoresearch.com) • Email: [info@zymoresearch.com](mailto:info@zymoresearch.com)



## ZYMO RESEARCH

*The Beauty of Science is to Make Things Simple*

Other kits will also transition using the below schedule:

<b>Product</b>	<b>Expected Date of Change</b>
<i>Quick-DNA Plant/Seed Miniprep Kit (Cat # D6020)</i>	<i>June 1<sup>st</sup>, 2025</i>
<i>Quick-RNA Plant Kit (Cat # R2024)</i>	<i>July 1<sup>st</sup>, 2025</i>
<i>Zymo Environ Water RNA Kit (Cat # R2042)</i>	<i>August 1<sup>st</sup>, 2025</i>
<i>Quick-RNA Fecal/Soil Microbe Microprep Kit (Cat # R2040)</i>	<i>September 1<sup>st</sup>, 2025</i>
<i>Quick-DNA Fecal/Soil Microbe Microprep Kit (Cat # D6012)</i>	<i>October 1<sup>st</sup>, 2025</i>
<i>Quick-DNA Fecal/Soil Microbe Miniprep Kit (Cat # D6010)</i>	<i>November 1<sup>st</sup>, 2025</i>
<i>Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat # D6110)</i>	<i>December 1<sup>st</sup>, 2025</i>

If you have further inquiries, please contact Zymo Research at [info@zymoresearch.com](mailto:info@zymoresearch.com) or [tech@zymoresearch.com](mailto:tech@zymoresearch.com).

We request that you share this notification with the relevant personnel in your organization.

Thank you for your continued partnership.

Brian D. Janssen, PhD  
International Distributor Account Manager  
Zymo Research Corporation