

Product Discontinuation Notification

Date: December 19, 2025

Product Change Notification: RD.PCN.139-A

Dear Valued Customer,

We are writing to inform you of an update to the previously announced transition timeline for the OneStep PCR Inhibitor Removal Kit.

Effective July 2025, Zymo Research began to decrease the supply of the OneStep PCR inhibitor removal kit (D6030) and Zymo-Spin III-HRC Filters (C1058-50) in an effort to fully transition to the new and improved format. We would like to notify you that the full transition to the updated format has been extended, and the previous format will now be available through June 30th, 2026.

The newly released and updated kit, [OneStep PCR Inhibitor Removal kit \(catalog number D6031\)](#), will now come with the new Zymo-Spin IV-IR HRC Filters (C1010-50), which serve the same function as the Zymo-Spin III-HRC Filters (C1058-50), but have a new design and chemical formulation, allowing it to remove PCR inhibitors more efficiently. In addition, the updated HRC IV-IR filters come in a wet slurry format, preventing column breakage and eliminating the need to add prep buffer prior to use.

D6030 OneStep PCR Inhibitor Removal kit

Zymo-Spin III-HRC Filters (C1058-50)



D6031 OneStep PCR Inhibitor Removal kit (NEW)

Zymo-Spin IV-IR HRC Filters (C1010-50)



Highlights

- Fast, one step procedure for cleaning impure samples prior to PCR, sequencing, NGS, reverse transcription (RT), etc.
- Efficiently eliminate polyphenolics from extracted DNA & RNA samples for PCR and short or long read sequencing.
- ≥ 80% recovery of PCR inhibitor-free DNA & RNA.
- Learn more about this new format: <https://www.zymoresearch.com/pages/onestep-pcr-inhibitor-removal-kit>

Orders for the previous format of the OneStep PCR inhibitor removal kit (D6030) and Zymo-Spin III-HRC Filters (C1058-50) will continue to be fulfilled until June 2026 or until inventories are depleted. If you have questions, need more information, or require assistance with transitioning, please do not hesitate to contact our technical support team at tech@zymoresearch.com.

Summary of change: Beginning July 2025, supply of the OneStep PCR Inhibitor Removal Kit (D6030) has been gradually reduced to support the transition to the updated OneStep PCR Inhibitor Removal Kit (D6031). The discontinuation date for the previous format has been extended from January 1st, 2026 to June 30th, 2026.

A copy of this notification should be immediately forwarded to the personnel in your organization affected by the change.

If you have any questions or concerns, please reach out to tech@zymoresearch.com.

Sincerely,

A handwritten signature in black ink, appearing to read 'Erica Carrasco', with a stylized, flowing script.

Erica Carrasco, MSc.
Scientific Collaboration Manager
Zymo Research Corporation