

Product Change Notification

Date: December 19, 2025

Product Change Notification #: RD.PCN.140-A

Dear Valued Distributor/Customer:

We are writing to inform you of an update to the previously announced transition timeline for the transition to the updated *OneStep* HRC Column.

Several of our environmental DNA and RNA extraction kits are in the process of transitioning from the Zymo-Spin III-HRC Filters (C1058-50) to the new updated version Zymo-Spin IV-IR HRC Filters (C1010-50). As part of our commitment to quality and continuous process improvements, we have integrated customer feedback and recommendations into the new design. The revised inhibitor removal columns feature a faster workflow, improved appearance, and more robust inhibitor removal performance.

What will change with this update:



- ✓ **Simplified Workflow and Improved Appearance:** HRC columns will now come with a slurry or “wet” format, deleting the need to add prep buffer prior to adding your sample.
- ✓ **Improved Performance:** Updated chemical formulation allows the HRC column to bind to more PCR inhibitors effectively cleaning up strongly inhibited samples in just one step!

The transition to the updated Zymo-Spin IV-IR HRC Filter format began with the Quick-RNA Plant Kit (Cat. No. R2024) on July 1, 2025. Customers who need to complete studies using the Zymo-Spin III-HRC Filter format with the Quick-RNA Plant Kit (R2024) were encouraged to place orders prior to that date.



Other kits will also transition using the below schedule:

Product	Expected Date of Change
Quick-DNA Plant/Seed Miniprep Kit (Cat # D6020)	June 1 st , 2025
Quick-RNA Plant Kit (Cat # R2024)	July 1 st , 2025
Zymo Environ Water RNA Kit (Cat # R2042)	August 1 st , 2025
Quick-RNA Fecal/Soil Microbe Microprep Kit (Cat # R2040)	September 1 st , 2025
Quick-DNA Fecal/Soil Microbe Microprep Kit (Cat # D6012)	October 1 st , 2025
Quick-DNA Fecal/Soil Microbe Miniprep Kit (Cat # D6010)	November 1 st , 2025
Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat # D6110)	Extended to June 30th, 2026

Please note that the transition timeline for Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat. No. D6110) has been extended, and the Zymo-Spin III-HRC Filter format will remain available through June 30th, 2026, or until inventory is depleted.

If you have further inquiries, please contact Zymo Research at info@zymoresearch.com or tech@zymoresearch.com.

We request that you share this notification with the relevant personnel in your organization.

Thank you for your continued partnership.

Erica Carrasco, Msc.
Scientific Collaboration Manager
Zymo Research Corp.