

## Product Change Notification

**Date:** April 23, 2025

**Product Change Notification #:** RD.PCN.140

Dear Valued Distributor/Customer:

This letter serves to inform you that several of our environmental DNA and RNA extraction kits will undergo a change of the Zymo-Spin III-HRC Filters (C1058-50) for the new updated version Zymo-Spin IV-IR HRC Filters (C1010-50). As part of our commitment to quality and continuous process improvements, we have integrated customer feedback and recommendations into the new design. The revised inhibitor removal columns feature a faster workflow, improved appearance, and more robust inhibitor removal performance.

What will change with this update:



- ✓ **Simplified Workflow and Improved Appearance:** HRC columns will now come with a slurry or “wet” format, deleting the need to add prep buffer prior to adding your sample.
- ✓ **Improved Performance:** Updated chemical formulation allows the HRC column to bind to more PCR inhibitors effectively cleaning up strongly inhibited samples in just one step!

The effort to transition to the new updated version continues with the kit **Quick-RNA Plant Kit (Cat # R2024)** starting **July 1<sup>st</sup>, 2025**. If you need to complete your study using *Quick-RNA Plant Kit* (Cat # R2024) using the Zymo-Spin III-HRC Filter format, please purchase prior to July 1st, 2025.



Other kits will also transition using the below schedule:

Product	Expected Date of Change
Quick-DNA Plant/Seed Miniprep Kit (Cat # D6020)	June 1 <sup>st</sup> , 2025
Quick-RNA Plant Kit (Cat # R2024)	July 1 <sup>st</sup> , 2025
Zymo Environ Water RNA Kit (Cat # R2042)	August 1 <sup>st</sup> , 2025
Quick-RNA Fecal/Soil Microbe Microprep Kit (Cat # R2040)	September 1 <sup>st</sup> , 2025
Quick-DNA Fecal/Soil Microbe Microprep Kit (Cat # D6012)	October 1 <sup>st</sup> , 2025
Quick-DNA Fecal/Soil Microbe Miniprep Kit (Cat # D6010)	November 1 <sup>st</sup> , 2025
Quick-DNA Fecal/Soil Microbe Midiprep Kit (Cat # D6110)	December 1 <sup>st</sup> , 2025

If you have further inquiries, please contact Zymo Research at [info@zymoresearch.com](mailto:info@zymoresearch.com) or [tech@zymoresearch.com](mailto:tech@zymoresearch.com).

We request that you share this notification with the relevant personnel in your organization.

Thank you for your continued partnership.

Brian D. Janssen, PhD  
International Distributor Account Manager  
Zymo Research Corporation