



Product Change Notification

Date: July 10th, 2026

Product Change Notification: RD.PCN.147

Dear Valued Distributor/Customer:

This letter serves to inform you of an important update to the protocol for the **Zymo-Seq SPLAT DNA Library Kit (D5464)**. These changes reflect our ongoing commitment to quality and ease of use. These changes are considered low impact and low risk, as they do not affect the performance of the products.

Below summarizes all updates included in this version release:

- **Reduced minimum input:** The kit now supports inputs as low as 1 ng (previously 10 ng), expanding its utility for low-input samples.
- **Added Library Bind Solution (D5464-A):** Improves the post-PCR clean-up efficiency. Post-ligation and post-PCR clean-up steps have been updated to reflect protocol changes and consolidated into their respective protocol sections to reflect this addition.

To facilitate a seamless transition:

- **Implementation Timeline:** This change is effective as of the date of this Product Change Notification. The updated kit configuration will be implemented beginning with the next product lot manufactured.
- **Remaining Inventory:** Existing inventory manufactured prior to the effective date should no longer be used or distributed, as it reflects the previous kit configuration and protocol, which have since been superseded. Please contact tech@zymoresearch.com with any questions regarding existing inventory.

If you have further inquiries, please contact Zymo Research at info@zymoresearch.com or tech@zymoresearch.com.

We request that you share this notification with the relevant personnel in your organization. Thank you for your continued partnership.

Erica Carrasco, MSc.
Scientific Collaboration Manager
Zymo Research Corp.